

TimeTagger

2.2.0

Generated by Doxygen 1.8.14

Contents

1	TimeTagger	1
2	Deprecated List	3
3	Module Index	5
3.1	Modules	5
4	Hierarchical Index	7
4.1	Class Hierarchy	7
5	Class Index	9
5.1	Class List	9
6	Module Documentation	11
6.1	TimeTagger backend	11
6.1.1	Detailed Description	11
6.2	base iterators	12
6.2.1	Detailed Description	13
6.2.2	Function Documentation	13
6.2.2.1	Dump()	13
6.2.2.2	TimeTagStream()	13

7 Class Documentation	15
7.1 <code>_Iterator</code> Class Reference	15
7.1.1 Detailed Description	17
7.1.2 Constructor & Destructor Documentation	17
7.1.2.1 <code>_Iterator()</code>	17
7.1.2.2 <code>~_Iterator()</code>	17
7.1.3 Member Function Documentation	17
7.1.3.1 <code>clear()</code>	17
7.1.3.2 <code>clear_impl()</code>	17
7.1.3.3 <code>getCaptureDuration()</code>	18
7.1.3.4 <code>getNewVirtualChannel()</code>	18
7.1.3.5 <code>isRunning()</code>	18
7.1.3.6 <code>lock()</code>	18
7.1.3.7 <code>next_impl()</code>	18
7.1.3.8 <code>registerChannel()</code>	19
7.1.3.9 <code>start()</code>	19
7.1.3.10 <code>startFor()</code>	19
7.1.3.11 <code>stop()</code>	20
7.1.3.12 <code>unlock()</code>	20
7.1.3.13 <code>unregisterChannel()</code>	20
7.1.4 Member Data Documentation	20
7.1.4.1 <code>channels_registered</code>	20
7.1.4.2 <code>running</code>	21
7.2 <code>Array< T ></code> Class Template Reference	21
7.3 Coincidence Class Reference	21
7.3.1 Detailed Description	22
7.3.2 Constructor & Destructor Documentation	22
7.3.2.1 <code>Coincidence()</code>	22
7.3.3 Member Function Documentation	22
7.3.3.1 <code>getChannel()</code>	22

7.4	Coincidences Class Reference	23
7.4.1	Detailed Description	23
7.4.2	Constructor & Destructor Documentation	24
7.4.2.1	Coincidences()	24
7.4.3	Member Function Documentation	24
7.4.3.1	getChannels()	24
7.4.3.2	next_impl()	24
7.5	Combiner Class Reference	25
7.5.1	Detailed Description	25
7.5.2	Constructor & Destructor Documentation	26
7.5.2.1	Combiner()	26
7.5.3	Member Function Documentation	26
7.5.3.1	clear_impl()	26
7.5.3.2	GET_DATA_1D()	26
7.5.3.3	getChannel()	27
7.5.3.4	next_impl()	27
7.6	Correlation Class Reference	27
7.6.1	Detailed Description	28
7.6.2	Constructor & Destructor Documentation	28
7.6.2.1	Correlation()	28
7.6.3	Member Function Documentation	29
7.6.3.1	clear_impl()	29
7.6.3.2	GET_DATA_1D() [1/2]	29
7.6.3.3	GET_DATA_1D() [2/2]	29
7.6.3.4	next_impl()	30
7.7	CountBetweenMarkers Class Reference	30
7.7.1	Detailed Description	31
7.7.2	Constructor & Destructor Documentation	31
7.7.2.1	CountBetweenMarkers()	31
7.7.3	Member Function Documentation	32

7.7.3.1	<code>clear_impl()</code>	32
7.7.3.2	<code>GET_DATA_1D()</code> [1/3]	32
7.7.3.3	<code>GET_DATA_1D()</code> [2/3]	32
7.7.3.4	<code>GET_DATA_1D()</code> [3/3]	33
7.7.3.5	<code>next_impl()</code>	33
7.7.3.6	<code>ready()</code>	33
7.8	Counter Class Reference	34
7.8.1	Detailed Description	34
7.8.2	Constructor & Destructor Documentation	34
7.8.2.1	<code>Counter()</code>	35
7.8.3	Member Function Documentation	36
7.8.3.1	<code>clear_impl()</code>	36
7.8.3.2	<code>GET_DATA_2D()</code>	36
7.8.3.3	<code>next_impl()</code>	36
7.8.3.4	<code>start()</code>	37
7.9	Countrate Class Reference	37
7.9.1	Detailed Description	38
7.9.2	Constructor & Destructor Documentation	38
7.9.2.1	<code>Countrate()</code>	38
7.9.3	Member Function Documentation	38
7.9.3.1	<code>clear_impl()</code>	38
7.9.3.2	<code>GET_DATA_1D()</code>	39
7.9.3.3	<code>next_impl()</code>	39
7.9.3.4	<code>start()</code>	39
7.10	DelayedChannel Class Reference	40
7.10.1	Detailed Description	40
7.10.2	Constructor & Destructor Documentation	40
7.10.2.1	<code>DelayedChannel()</code>	40
7.10.3	Member Function Documentation	41
7.10.3.1	<code>getChannel()</code>	41

7.10.3.2	next_impl()	41
7.10.3.3	setDelay()	42
7.10.3.4	start()	42
7.11	Dump Class Reference	42
7.11.1	Detailed Description	43
7.11.2	Constructor & Destructor Documentation	43
7.11.2.1	~Dump()	43
7.11.3	Member Function Documentation	43
7.11.3.1	clear_impl()	43
7.11.3.2	next_impl()	43
7.12	Event Struct Reference	44
7.13	Flim Class Reference	44
7.13.1	Detailed Description	45
7.13.2	Constructor & Destructor Documentation	45
7.13.2.1	Flim()	45
7.14	FrequencyMultiplier Class Reference	45
7.14.1	Detailed Description	46
7.14.2	Constructor & Destructor Documentation	46
7.14.2.1	FrequencyMultiplier()	46
7.14.3	Member Function Documentation	47
7.14.3.1	next_impl()	47
7.15	GatedChannel Class Reference	47
7.15.1	Detailed Description	48
7.15.2	Constructor & Destructor Documentation	48
7.15.2.1	GatedChannel()	48
7.15.3	Member Function Documentation	49
7.15.3.1	getChannel()	49
7.15.3.2	next_impl()	49
7.16	Histogram Class Reference	49
7.16.1	Detailed Description	50

7.16.2	Constructor & Destructor Documentation	50
7.16.2.1	Histogram()	50
7.17	HistogramLogBins Class Reference	51
7.17.1	Detailed Description	51
7.17.2	Constructor & Destructor Documentation	51
7.17.2.1	HistogramLogBins()	52
7.17.3	Member Function Documentation	52
7.17.3.1	clear_impl()	52
7.17.3.2	next_impl()	52
7.18	Iterator Class Reference	53
7.18.1	Detailed Description	54
7.18.2	Constructor & Destructor Documentation	54
7.18.2.1	Iterator()	54
7.18.3	Member Function Documentation	54
7.18.3.1	clear_impl()	54
7.18.3.2	next()	54
7.18.3.3	next_impl()	55
7.18.3.4	size()	55
7.19	LED Union Reference	55
7.20	PRBS Class Reference	56
7.21	Scope Class Reference	56
7.21.1	Constructor & Destructor Documentation	57
7.21.1.1	Scope()	57
7.21.2	Member Function Documentation	57
7.21.2.1	clear_impl()	57
7.21.2.2	next_impl()	58
7.21.2.3	start()	58
7.22	StartStop Class Reference	58
7.22.1	Detailed Description	59
7.22.2	Constructor & Destructor Documentation	59

7.22.2.1	StartStop()	59
7.22.3	Member Function Documentation	60
7.22.3.1	clear_impl()	60
7.22.3.2	next_impl()	60
7.22.3.3	start()	60
7.23	SynchronizedMeasurments Class Reference	61
7.23.1	Detailed Description	61
7.23.2	Constructor & Destructor Documentation	61
7.23.2.1	SynchronizedMeasurments()	61
7.23.3	Member Function Documentation	62
7.23.3.1	clear()	62
7.23.3.2	registerMeasurement()	62
7.23.3.3	start()	62
7.23.3.4	stop()	62
7.24	Tag Struct Reference	63
7.24.1	Detailed Description	63
7.24.2	Member Data Documentation	63
7.24.2.1	channel	63
7.24.2.2	overflow	63
7.24.2.3	time	64
7.25	TimeDifferences Class Reference	64
7.25.1	Detailed Description	65
7.25.2	Constructor & Destructor Documentation	65
7.25.2.1	TimeDifferences()	65
7.25.3	Member Function Documentation	65
7.25.3.1	clear_impl()	66
7.25.3.2	GET_DATA_2D()	66
7.25.3.3	getCounts()	66
7.25.3.4	next_impl()	66
7.25.3.5	ready()	67

7.25.3.6	setMaxCounts()	67
7.25.3.7	start()	67
7.26	TimeTagger Class Reference	68
7.26.1	Detailed Description	70
7.26.2	Constructor & Destructor Documentation	70
7.26.2.1	TimeTagger()	70
7.26.2.2	~TimeTagger()	70
7.26.3	Member Function Documentation	70
7.26.3.1	autoCalibration()	70
7.26.3.2	clearOverflows()	71
7.26.3.3	GET_DATA_2D() [1/2]	71
7.26.3.4	GET_DATA_2D() [2/2]	71
7.26.3.5	getChannelNumberScheme()	71
7.26.3.6	getChannels()	72
7.26.3.7	getConditionalFilterFiltered()	72
7.26.3.8	getConditionalFilterTrigger()	72
7.26.3.9	getDACRange()	72
7.26.3.10	getDeadtime()	72
7.26.3.11	getEventDivider()	73
7.26.3.12	getFilter()	73
7.26.3.13	getHardwareBufferSize()	73
7.26.3.14	getInputDelay()	73
7.26.3.15	getInvertedChannel()	74
7.26.3.16	getModel()	74
7.26.3.17	getNormalization()	74
7.26.3.18	getOverflows()	74
7.26.3.19	getOverflowsAndClear()	74
7.26.3.20	getPcbVersion()	75
7.26.3.21	getPsPerClock()	75
7.26.3.22	getSensorData()	75

7.26.3.23	getSerial()	75
7.26.3.24	getTestSignal()	75
7.26.3.25	getTriggerLevel()	76
7.26.3.26	isUnusedChannel()	76
7.26.3.27	registerChannel()	76
7.26.3.28	reset()	76
7.26.3.29	runSynchronized()	77
7.26.3.30	setConditionalFilter()	77
7.26.3.31	setDeadtime()	77
7.26.3.32	setEventDivider()	78
7.26.3.33	setFilter()	78
7.26.3.34	setHardwareBufferSize()	79
7.26.3.35	setInputDelay()	79
7.26.3.36	setLED()	79
7.26.3.37	setNormalization()	79
7.26.3.38	setTestSignal()	80
7.26.3.39	setTestSignalDivider()	80
7.26.3.40	setTriggerLevel()	80
7.26.3.41	sync()	81
7.26.3.42	unregisterChannel()	81
7.27	TimetaggerFPGA Class Reference	81
7.28	TimeTaggerModel Class Reference	82
7.29	TimeTagStream Class Reference	82
7.29.1	Detailed Description	83
7.29.2	Constructor & Destructor Documentation	83
7.29.2.1	~TimeTagStream()	83
7.29.3	Member Function Documentation	83
7.29.3.1	clear_impl()	83
7.29.3.2	getCounts()	84
7.29.3.3	next_impl()	84
7.30	TimeTagStreamBuffer Class Reference	84

Chapter 1

TimeTagger

backend for [TimeTagger](#), an OpalKelly based single photon counting library

Author

Markus Wick wickmarkus@web.de

Helmut Fedder helmut.fedder@gmail.com

[TimeTagger](#) provides an easy to use and cost effective hardware solution for time-resolved single photon counting applications.

This document describes the C++ native interface to the [TimeTagger](#) device.

Chapter 2

Deprecated List

Member `TimeTagger::getChannels` ()=0

Use `getChannelList` instead.

Member `TimeTagger::getFilter` ()=0

use `getConditionalFilter*`

Member `TimeTagger::setFilter` (bool state)=0

use `setConditionalFilter`

Chapter 3

Module Index

3.1 Modules

Here is a list of all modules:

TimeTagger backend	11
base iterators	12

Chapter 4

Hierarchical Index

4.1 Class Hierarchy

This inheritance list is sorted roughly, but not completely, alphabetically:

_Iterator	15
Coincidences	23
Coincidence	21
Combiner	25
Correlation	27
CountBetweenMarkers	30
Counter	34
Countrate	37
DelayedChannel	40
Dump	42
FrequencyMultiplier	45
GatedChannel	47
HistogramLogBins	51
Iterator	53
Scope	56
StartStop	58
TimeDifferences	64
TimeTagStream	82
Array< T >	21
Event	44
Flim	44
Histogram	49
LED	55
PRBS	56
SynchronizedMeasurments	61
Tag	63
TimeTagger	68
TimetaggerFPGA	81
TimeTaggerModel	82
TimeTagStreamBuffer	84

Chapter 5

Class Index

5.1 Class List

Here are the classes, structs, unions and interfaces with brief descriptions:

_Iterator	Base class for all iterators	15
Array< T >	21
Coincidence	Coincidence monitor for one or more channel groups	21
Coincidences	Coincidence monitor for one or more channel groups	23
Combiner	Combine some channels in a virtual channel which has a tick for each tick in the input channels	25
Correlation	Cross-correlation between two channels	27
CountBetweenMarkers	Simple counter where external marker signals determine the bins	30
Counter	Simple counter on one or more channels	34
Countrate	Count rate on one or more channels	37
DelayedChannel	Simple delayed queue	40
Dump	Dump all time tags to a file	42
Event	44
Flim	Fluorescence lifetime imaging	44
FrequencyMultiplier	The signal of an input channel is scaled up to a higher frequency according to the multiplier passed as a parameter	45
GatedChannel	An input channel is gated by a gate channel	47
Histogram	Accumulate time differences into a histogram	49
HistogramLogBins	Accumulate time differences into a histogram with logarithmic increasing bin sizes	51
Iterator	Simple event queue	53

LED	55
PRBS	56
Scope	56
StartStop	
Simple start-stop measurement	58
SynchronizedMeasurements	
Start, stop and clear several measurements synchronized	61
Tag	
Single event on a channel	63
TimeDifferences	
Accumulates the time differences between clicks on two channels in one or more histograms	64
TimeTagger	
Backend for the TimeTagger	68
TimeTaggerFPGA	81
TimeTaggerModel	82
TimeTagStream	
Access the time tag stream	82
TimeTagStreamBuffer	84

Chapter 6

Module Documentation

6.1 TimeTagger backend

the timetagger kernel

Classes

- class [TimeTagger](#)
backend for the [TimeTagger](#).
- struct [Tag](#)
a single event on a channel
- class [_Iterator](#)
Base class for all iterators.

6.1.1 Detailed Description

the timetagger kernel

6.2 base iterators

base iterators for photon counting applications

Classes

- class [Combiner](#)
Combine some channels in a virtual channel which has a tick for each tick in the input channels.
- class [CountBetweenMarkers](#)
a simple counter where external marker signals determine the bins
- class [Counter](#)
a simple counter on one or more channels
- class [Coincidences](#)
a coincidence monitor for one or more channel groups
- class [Coincidence](#)
a coincidence monitor for one or more channel groups
- class [Countrate](#)
count rate on one or more channels
- class [DelayedChannel](#)
a simple delayed queue
- class [GatedChannel](#)
An input channel is gated by a gate channel.
- class [FrequencyMultiplier](#)
The signal of an input channel is scaled up to a higher frequency according to the multiplier passed as a parameter.
- class [Iterator](#)
a simple event queue
- class [StartStop](#)
simple start-stop measurement
- class [TimeDifferences](#)
Accumulates the time differences between clicks on two channels in one or more histograms.
- class [Histogram](#)
Accumulate time differences into a histogram.
- class [HistogramLogBins](#)
Accumulate time differences into a histogram with logarithmic increasing bin sizes.
- class [Flim](#)
Fluorescence lifetime imaging.
- class [Correlation](#)
cross-correlation between two channels
- class [Scope](#)

Functions

- [TimeTagStream::TimeTagStream](#) ([TimeTagger](#) *tagger, int n_max_events, const std::vector< channel_t > &channels=std::vector< channel_t >())
constructor of a [TimeTagStream](#) thread
- [Dump::Dump](#) ([TimeTagger](#) *tagger, std::string filename, size_t max_tags, const std::vector< channel_t > &channels=std::vector< channel_t >())
constructor of a [Dump](#) thread

6.2.1 Detailed Description

base iterators for photon counting applications

6.2.2 Function Documentation

6.2.2.1 Dump()

```
Dump::Dump (
    TimeTagger * tagger,
    std::string filename,
    size_t max_tags,
    const std::vector< channel_t > & channels = std::vector<channel_t>() )
```

constructor of a [Dump](#) thread

Parameters

<i>tagger</i>	reference to a TimeTagger
<i>filename</i>	name of the file to dump to
<i>max_tags</i>	stop after this number of tags has been dumped
<i>channels</i>	channels which are dumped to the file (when empty or not passed all active channels are dumped)

6.2.2.2 TimeTagStream()

```
TimeTagStream::TimeTagStream (
    TimeTagger * tagger,
    int n_max_events,
    const std::vector< channel_t > & channels = std::vector<channel_t>() )
```

constructor of a [TimeTagStream](#) thread

Gives access to the time tag stream

Parameters

<i>tagger</i>	reference to a TimeTagger
<i>n_max_events</i>	maximum number of tags stored
<i>channels</i>	channels which are dumped to the file (when empty or not passed all active channels are dumped)

Chapter 7


Class Documentation

7.1 `_Iterator` Class Reference

Base class for all iterators.

```
#include <TimeTagger.h>
```

Inheritance diagram for `_Iterator`:



Public Member Functions

- virtual `~_Iterator` ()
destructor
- virtual void `start` ()
start the iterator
- void `startFor` (timestamp_t capture_duration, bool `clear`=true)
start the iterator, and stops it after the capture_duration
- virtual void `stop` ()
stop the iterator
- void `clear` ()
clear `Iterator` state.
- bool `isRunning` ()
query the `Iterator` state.
- timestamp_t `getCaptureDuration` ()
query the evaluation time

Protected Member Functions

- `_Iterator` (TimeTagger *tagger)
standard constructor
- void `registerChannel` (channel_t channel)
register a channel
- void `unregisterChannel` (channel_t channel)
unregister a channel
- channel_t `getNewVirtualChannel` ()
allocate a new virtual output channel for this iterator
- virtual void `clear_impl` ()
clear `Iterator` state.
- void `lock` ()
acquire update lock
- void `unlock` ()
release update lock
- virtual bool `next_impl` (std::vector< Tag > &incoming_tags, timestamp_t begin_time, timestamp_t end_time)=0
update iterator state

Protected Attributes

- std::set< channel_t > `channels_registered`
list of channels used by the iterator
- bool `running`
running state of the iterator
- TimeTagger * `tagger`
- timestamp_t `capture_duration`

Friends

- class `_TimeTagger`

7.1.1 Detailed Description

Base class for all iterators.

7.1.2 Constructor & Destructor Documentation

7.1.2.1 [_Iterator\(\)](#)

```
_Iterator::_Iterator (
    TimeTagger * tagger ) [protected]
```

standard constructor

will register with the [TimeTagger](#) backend.

7.1.2.2 [~_Iterator\(\)](#)

```
virtual _Iterator::~~_Iterator ( ) [virtual]
```

destructor

will stop and unregister prior finalization.

7.1.3 Member Function Documentation

7.1.3.1 [clear\(\)](#)

```
void _Iterator::clear ( )
```

clear [iterator](#) state.

7.1.3.2 [clear_impl\(\)](#)

```
virtual void _Iterator::clear_impl ( ) [inline], [protected], [virtual]
```

clear [iterator](#) state.

Each [iterator](#) should implement the [clear_impl\(\)](#) method to reset its internal state. The [clear_impl\(\)](#) function is guarded by the update lock.

Reimplemented in [Scope](#), [Correlation](#), [HistogramLogBins](#), [TimeDifferences](#), [StartStop](#), [Dump](#), [TimeTagStream](#), [Iterator](#), [Countrate](#), [Counter](#), [CountBetweenMarkers](#), and [Combiner](#).

7.1.3.3 `getCaptureDuration()`

```
timestamp_t _Iterator::getCaptureDuration ( )
```

query the evaluation time

Query the total capture duration since the last call to `.`. This might have a wrong amount of time if there were some overflows within this range.

Returns

capture duration of the data

7.1.3.4 `getNewVirtualChannel()`

```
channel_t _Iterator::getNewVirtualChannel ( ) [protected]
```

allocate a new virtual output channel for this iterator

7.1.3.5 `isRunning()`

```
bool _Iterator::isRunning ( )
```

query the [Iterator](#) state.

Fetches if this iterator is running.

7.1.3.6 `lock()`

```
void _Iterator::lock ( ) [protected]
```

acquire update lock

All mutable operations on a iterator are guarded with an update mutex. Implementers are adviced to [lock\(\)](#) an iterator, whenever internal state is queried or changed.

7.1.3.7 `next_impl()`

```
virtual bool _Iterator::next_impl (
    std::vector< Tag > & incoming_tags,
    timestamp_t begin_time,
    timestamp_t end_time ) [protected], [pure virtual]
```

update iterator state

Each [Iterator](#) must implement the [next_impl\(\)](#) method. The [next_impl\(\)](#) function is guarded by the update lock.

The backend delivers each [Tag](#) on each registered channel to this callback function.

Parameters

<i>list</i>	block of events
<i>start_time</i>	earliest event in the block
<i>end_time</i>	start_time of the next block, not including in this block

Returns

if the content of this block was modified

Implemented in [Scope](#), [Correlation](#), [HistogramLogBins](#), [TimeDifferences](#), [StartStop](#), [Dump](#), [TimeTagStream](#), [Iterator](#), [FrequencyMultiplier](#), [GatedChannel](#), [DelayedChannel](#), [Countrate](#), [Coincidences](#), [Counter](#), [CountBetweenMarkers](#), and [Combiner](#).

7.1.3.8 registerChannel()

```
void _Iterator::registerChannel (
    channel_t channel ) [protected]
```

register a channel

Only channels registered by any iterator attached to a backend are delivered over the usb.

Parameters

<i>channel</i>	the channel
----------------	-------------

7.1.3.9 start()

```
virtual void _Iterator::start ( ) [virtual]
```

start the iterator

The default behavior for iterators is to start automatically on creation.

Reimplemented in [Scope](#), [TimeDifferences](#), [StartStop](#), [DelayedChannel](#), [Countrate](#), and [Counter](#).

7.1.3.10 startFor()

```
void _Iterator::startFor (
    timestamp_t capture_duration,
    bool clear = true )
```

start the iterator, and stops it after the capture_duration

Parameters

<i>capture_duration</i>	capture duration until the measurement is stopped
<i>clear</i>	resets the data acquired

When the startFor is called before the previous measurement has ended and the clear parameter is set to false, then the passed capture_duration will be added on top to the current max_capture_duration

7.1.3.11 stop()

```
virtual void _Iterator::stop ( ) [virtual]
```

stop the iterator

The iterator is put into the STOPPED state, but will still be registered with the backend.

7.1.3.12 unlock()

```
void _Iterator::unlock ( ) [protected]
```

release update lock

see [lock\(\)](#)

7.1.3.13 unregisterChannel()

```
void _Iterator::unregisterChannel (
    channel_t channel ) [protected]
```

unregister a channel

Parameters

<i>channel</i>	the channel
----------------	-------------

7.1.4 Member Data Documentation**7.1.4.1 channels_registered**

```
std::set<channel_t> _Iterator::channels_registered [protected]
```

list of channels used by the iterator

7.1.4.2 running

```
bool _Iterator::running [protected]
```

running state of the iterator

The documentation for this class was generated from the following file:

- TimeTagger.h

7.2 Array< T > Class Template Reference

Public Member Functions

- **Array** (std::vector< int > dims)
- int **getDims** ()
- int **getSize** (int dim)
- size_t **getLength** ()
- void **get** (T *fixed_output)
- T * **get** ()
- **Array** (std::vector< int > dims)
- int **getDims** ()
- int **getSize** (int dim)
- size_t **getLength** ()
- void **get** (T *fixed_output)
- T * **get** ()

The documentation for this class was generated from the following files:


- CSharpWrapper.h
- JavaWrapper.h

7.3 Coincidence Class Reference

a coincidence monitor for one or more channel groups

```
#include <Iterators.h>
```

Inheritance diagram for Coincidence:



Public Member Functions

- [Coincidence](#) ([TimeTagger](#) *tagger, std::vector< channel_t > channels, timestamp_t coincidenceWindow=1000)
construct a coincidence
- channel_t [getChannel](#) ()
virtual channel which contains the coincidences

Additional Inherited Members

7.3.1 Detailed Description

a coincidence monitor for one or more channel groups

Monitor coincidences for a given channel groups passed by the constructor. A coincidence is event is detected when all slected channels have a click within the given coincidenceWindow [ps] The coincidence will create a virtual events on a virtual channel with the channel number provided by [getChannel\(\)](#). For multiple coincidence channel combinations use the class [Coincidences](#) which outperforms multiple instances of Conincidence.

7.3.2 Constructor & Destructor Documentation

7.3.2.1 Coincidence()

```
Coincidence::Coincidence (
    TimeTagger * tagger,
    std::vector< channel_t > channels,
    timestamp_t coincidenceWindow = 1000 ) [inline]
```

construct a coincidence

Parameters

<i>tagger</i>	reference to a TimeTagger
<i>channels</i>	vector of channels to match
<i>window</i>	max distance between all clicks for a coincidence [ps]

7.3.3 Member Function Documentation

7.3.3.1 getChannel()

```
channel_t Coincidence::getChannel ( ) [inline]
```

virtual channel which contains the coincidences

The documentation for this class was generated from the following file:


- `Iterators.h`

7.4 Coincidences Class Reference

a coincidence monitor for one or more channel groups

```
#include <Iterators.h>
```

Inheritance diagram for Coincidences:



Public Member Functions

- `Coincidences` (`TimeTagger` *tagger, `std::vector< std::vector< channel_t >>` coincidenceGroups, `timestamp_t` coincidenceWindow)
construct a `Coincidences`
- `std::vector< channel_t >` `getChannels` ()
fetches the block of virtual channels for those coincidence groups
- `void` `setCoincidenceWindow` (`timestamp_t` coincidenceWindow)

Protected Member Functions

- `bool` `next_impl` (`std::vector< Tag >` &incoming_tags, `timestamp_t` begin_time, `timestamp_t` end_time) override
update iterator state

Additional Inherited Members

7.4.1 Detailed Description

a coincidence monitor for one or more channel groups

Monitor coincidences for given coincidence groups passed by the constructor. A coincidence is hereby defined as for a given coincidence group a) the incoming is part of this group b) at least tag arrived within the coincidence↔ Window [ps] for all other channels of this coincidence group Each coincidence will create a virtual event. The block of event IDs for those coincidence group can be fetched.

7.4.2 Constructor & Destructor Documentation

7.4.2.1 Coincidences()

```
Coincidences::Coincidences (
    TimeTagger * tagger,
    std::vector< std::vector< channel_t >> coincidenceGroups,
    timestamp_t coincidenceWindow )
```

construct a [Coincidences](#)

Parameters

<i>tagger</i>	reference to a TimeTagger
<i>coincidenceGroups</i>	a vector of channels defining the coincidences
<i>coincidenceWindow</i>	the size of the coincidence window in picoseconds

7.4.3 Member Function Documentation

7.4.3.1 getChannels()

```
std::vector< channel_t > Coincidences::getChannels ( )
```

fetches the block of virtual channels for those coincidence groups

7.4.3.2 next_impl()

```
bool Coincidences::next_impl (
    std::vector< Tag > & incoming_tags,
    timestamp_t begin_time,
    timestamp_t end_time ) [override], [protected], [virtual]
```

update iterator state

Each [Iterator](#) must implement the [next_impl\(\)](#) method. The [next_impl\(\)](#) function is guarded by the update lock.

The backend delivers each [Tag](#) on each registered channel to this callback function.

Parameters

<i>list</i>	block of events
<i>start_time</i>	earliest event in the block
<i>end_time</i>	start_time of the next block, not including in this block

Returns

if the content of this block was modified

Implements [_Iterator](#).

The documentation for this class was generated from the following files:


- Iterators.h
- Coincidence.cpp

7.5 Combiner Class Reference

Combine some channels in a virtual channel which has a tick for each tick in the input channels.

```
#include <Iterators.h>
```

Inheritance diagram for Combiner:

**Public Member Functions**

- [Combiner](#) ([TimeTagger](#) *tagger, std::vector< [channel_t](#) > channels)
construct a combiner
- [GET_DATA_1D](#) (getData, long long, array_out,)
get sum of counts
- [channel_t](#) [getChannel](#) ()
the new virtual channel

Protected Member Functions

- bool [next_impl](#) (std::vector< [Tag](#) > &incoming_tags, timestamp_t begin_time, timestamp_t end_time) override
update iterator state
- void [clear_impl](#) () override
clear [Iterator](#) state.

Additional Inherited Members

7.5.1 Detailed Description

Combine some channels in a virtual channel which has a tick for each tick in the input channels.

This iterator can be used to get aggregation channels, eg if you want to monitor the countrate of the sum of two channels.

7.5.2 Constructor & Destructor Documentation

7.5.2.1 Combiner()

```
Combiner::Combiner (
    TimeTagger * tagger,
    std::vector< channel_t > channels )
```

construct a combiner

Parameters

<i>tagger</i>	reference to a TimeTagger
<i>channels</i>	vector of channels to combine

7.5.3 Member Function Documentation

7.5.3.1 clear_impl()

```
void Combiner::clear_impl ( ) [override], [protected], [virtual]
```

clear [Iterator](#) state.

Each [Iterator](#) should implement the [clear_impl\(\)](#) method to reset its internal state. The [clear_impl\(\)](#) function is guarded by the update lock.

Reimplemented from [_Iterator](#).

7.5.3.2 GET_DATA_1D()

```
Combiner::GET_DATA_1D (
    getData ,
    long long,
    array_out )
```

get sum of counts

For reference, this iterators sums up how much ticks are generated because of which input channel. So this functions returns an array with one value per input channel.

7.5.3.3 `getChannel()`

```
channel_t Combiner::getChannel ( )
```

the new virtual channel

This function returns the new allocated virtual channel. It can be used now in any new iterator.

7.5.3.4 `next_impl()`

```
bool Combiner::next_impl (
    std::vector< Tag > & incoming_tags,
    timestamp_t begin_time,
    timestamp_t end_time ) [override], [protected], [virtual]
```

update iterator state

Each `Iterator` must implement the `next_impl()` method. The `next_impl()` function is guarded by the update lock.

The backend delivers each `Tag` on each registered channel to this callback function.

Parameters

<i>list</i>	block of events
<i>start_time</i>	earliest event in the block
<i>end_time</i>	start_time of the next block, not including in this block

Returns

if the content of this block was modified

Implements `_Iterator`.

The documentation for this class was generated from the following files:


- Iterators.h
- Combiner.cpp

7.6 Correlation Class Reference

cross-correlation between two channels

```
#include <Iterators.h>
```

Inheritance diagram for `Correlation`:



Public Member Functions

- [Correlation](#) ([TimeTagger](#) *tagger, channel_t channel_1, channel_t channel_2=CHANNEL_UNUSED, timestamp_t binwidth=1000, int n_bins=1000)
constructor of a correlation measurement
- [GET_DATA_1D](#) (getData, int, array_out,)
get result data
- [GET_DATA_1D](#) (getDataNormalized, double, array_out,)
get result data - normalized such that a perfectly uncorrelated signals would be flat at a height of one
- [GET_DATA_1D](#) (getIndex, timestamp_t, array_out,)

Protected Member Functions

- bool [next_impl](#) (std::vector< [Tag](#) > &incoming_tags, timestamp_t begin_time, timestamp_t end_time) override
update iterator state
- void [clear_impl](#) () override
clear [Iterator](#) state.

Additional Inherited Members

7.6.1 Detailed Description

cross-correlation between two channels

Accumulates time differences between clicks on two channels into a histogram, where all ticks are considered both as start and stop clicks and both positive and negative time differences are considered. The histogram is determined by the number of bins and the binwidth, which are used both for the positive and the negative histogram range (i.e., length of the histogram is $2*n_bins+1$).

7.6.2 Constructor & Destructor Documentation

7.6.2.1 Correlation()

```
Correlation::Correlation (
    TimeTagger * tagger,
    channel_t channel_1,
    channel_t channel_2 = CHANNEL_UNUSED,
    timestamp_t binwidth = 1000,
    int n_bins = 1000 )
```

constructor of a correlation measurement

Parameters

<i>tagger</i>	reference to a TimeTagger
<i>channel_1</i>	first click channel
<i>channel_2</i>	second click channel
<i>binwidth</i>	width of one histogram bin in ps
<i>n_bins</i>	the number of bins in the resulting histogram is $2*n_bins+1$

7.6.3 Member Function Documentation

7.6.3.1 clear_impl()

```
void Correlation::clear_impl ( ) [override], [protected], [virtual]
```

clear [Iterator](#) state.

Each [Iterator](#) should implement the [clear_impl\(\)](#) method to reset its internal state. The [clear_impl\(\)](#) function is guarded by the update lock.

Reimplemented from [_Iterator](#).

7.6.3.2 GET_DATA_1D() [1/2]

```
Correlation::GET_DATA_1D (
    getData ,
    int ,
    array_out )
```

get result data

Parameters

<i>ARGOUTVIEWM_ARRAY2</i>	pointer receiving pointer to data
<i>DIM1</i>	pointer receiving first dimension
<i>DIM2</i>	pointer receiving second dimension

7.6.3.3 GET_DATA_1D() [2/2]

```
Correlation::GET_DATA_1D (
    getDataNormalized ,
    double ,
    array_out )
```

get result data - normalized such that a perfectly uncorrelated signals would be flat at a height of one

Parameters

<i>ARGOUTVIEWM_ARRAY2</i>	pointer receiving pointer to data
<i>DIM1</i>	pointer receiving first dimension
<i>DIM2</i>	pointer receiving second dimension

7.6.3.4 next_impl()

```
bool Correlation::next_impl (
    std::vector< Tag > & incoming_tags,
    timestamp_t begin_time,
    timestamp_t end_time ) [override], [protected], [virtual]
```

update iterator state

Each [Iterator](#) must implement the [next_impl\(\)](#) method. The [next_impl\(\)](#) function is guarded by the update lock.

The backend delivers each [Tag](#) on each registered channel to this callback function.

Parameters

<i>list</i>	block of events
<i>start_time</i>	earliest event in the block
<i>end_time</i>	start_time of the next block, not including in this block

Returns

if the content of this block was modified

Implements [_Iterator](#).

The documentation for this class was generated from the following files:


- Iterators.h
- TimeDifferences.cpp

7.7 CountBetweenMarkers Class Reference

a simple counter where external marker signals determine the bins

```
#include <Iterators.h>
```

Inheritance diagram for CountBetweenMarkers:



Public Member Functions

- [CountBetweenMarkers](#) ([TimeTagger](#) *tagger, channel_t click_channel, channel_t begin_channel, channel_t end_channel=CHANNEL_UNUSED, int n_values=1000)
constructor of [CountBetweenMarkers](#)
- bool [ready](#) ()
tbd
- [GET_DATA_1D](#) (getData, int, array_out,)
tbd
- [GET_DATA_1D](#) (getBinWidths, timestamp_t, array_out,)
fetches the widths of each bins
- [GET_DATA_1D](#) (getIndex, timestamp_t, array_out,)
fetches the starting time of each bin

Protected Member Functions

- bool [next_impl](#) (std::vector< [Tag](#) > &incoming_tags, timestamp_t begin_time, timestamp_t end_time) override
update iterator state
- void [clear_impl](#) () override
clear [Iterator](#) state.

Additional Inherited Members

7.7.1 Detailed Description

a simple counter where external marker signals determine the bins

[Counter](#) with external signals that trigger beginning and end of each counter accumulation. This can be used to implement counting triggered by a pixel clock and gated counting. The thread waits for the first time tag on the 'begin_channel', then begins counting time tags on the 'click_channel'. It ends counting when a tag on the 'end_channel' is detected.

7.7.2 Constructor & Destructor Documentation

7.7.2.1 CountBetweenMarkers()

```
CountBetweenMarkers::CountBetweenMarkers (
    TimeTagger * tagger,
    channel_t click_channel,
    channel_t begin_channel,
    channel_t end_channel = CHANNEL_UNUSED,
    int n_values = 1000 )
```

constructor of [CountBetweenMarkers](#)

Parameters

<i>tagger</i>	reference to a TimeTagger
<i>click_channel</i>	channel that increases the count
<i>begin_channel</i>	channel that triggers beginning of counting and stepping to the next value
<i>end_channel</i>	channel that triggers end of counting
<i>n_values</i>	the number of counter values to be stored

7.7.3 Member Function Documentation

7.7.3.1 `clear_impl()`

```
void CountBetweenMarkers::clear_impl ( ) [override], [protected], [virtual]
```

clear [Iterator](#) state.

Each [Iterator](#) should implement the `clear_impl()` method to reset its internal state. The `clear_impl()` function is guarded by the update lock.

Reimplemented from [_Iterator](#).

7.7.3.2 `GET_DATA_1D()` [1/3]

```
CountBetweenMarkers::GET_DATA_1D (
    getData ,
    int ,
    array_out )
```

tbd

Parameters

<i>ARGOUTVIEWM_ARRAY1</i>	
<i>DIM1</i>	

7.7.3.3 `GET_DATA_1D()` [2/3]

```
CountBetweenMarkers::GET_DATA_1D (
    getBinWidths ,
    timestamp_t ,
    array_out )
```

fetches the widths of each bins

7.7.3.4 GET_DATA_1D() [3/3]

```
CountBetweenMarkers::GET_DATA_1D (
    getIndex ,
    timestamp_t ,
    array_out )
```

fetches the starting time of each bin

7.7.3.5 next_impl()

```
bool CountBetweenMarkers::next_impl (
    std::vector< Tag > & incoming_tags,
    timestamp_t begin_time,
    timestamp_t end_time ) [override], [protected], [virtual]
```

update iterator state

Each [Iterator](#) must implement the [next_impl\(\)](#) method. The [next_impl\(\)](#) function is guarded by the update lock.

The backend delivers each [Tag](#) on each registered channel to this callback function.

Parameters

<i>list</i>	block of events
<i>start_time</i>	earliest event in the block
<i>end_time</i>	start_time of the next block, not including in this block

Returns

if the content of this block was modified

Implements [_Iterator](#).

7.7.3.6 ready()

```
bool CountBetweenMarkers::ready ( )
```

tbd

The documentation for this class was generated from the following files:


- Iterators.h
- CountBetweenMarkers.cpp

7.8 Counter Class Reference

a simple counter on one or more channels

```
#include <Iterators.h>
```

Inheritance diagram for Counter:



Public Member Functions

- `Counter` (`TimeTagger` *tagger, `std::vector< channel_t >` channels, `timestamp_t` binwidth=1000000000, `int` n_values=1)
construct a counter
- `void start` () override
start the iterator
- `GET_DATA_2D` (`getData`, `int`, `array_out`,)
get counts
- `GET_DATA_1D` (`getIndex`, `timestamp_t`, `array_out`,)

Protected Member Functions

- `bool next_impl` (`std::vector< Tag >` &incoming_tags, `timestamp_t` begin_time, `timestamp_t` end_time) override
update iterator state
- `void clear_impl` () override
clear `Iterator` state.

Additional Inherited Members

7.8.1 Detailed Description

a simple counter on one or more channels

`Counter` with fixed binwidth and circular buffer output. This class is suitable to generate a time trace of the count rate on one or more channels. The thread repeatedly counts clicks on a single channel over a given time interval and stores the results in a two-dimensional array. The array is treated as a circular buffer. I.e., once the array is full, each new value shifts all previous values one element to the left.

7.8.2 Constructor & Destructor Documentation

7.8.2.1 Counter()

```
Counter::Counter (
    TimeTagger * tagger,
    std::vector< channel_t > channels,
    timestamp_t binwidth = 1000000000,
    int n_values = 1 )
```

construct a counter

Parameters

<i>tagger</i>	reference to a TimeTagger
<i>channels</i>	channels to count on
<i>binwidth</i>	counts are accumulated for binwidth picoseconds
<i>n_values</i>	number of counter values stored (for each channel)

7.8.3 Member Function Documentation

7.8.3.1 clear_impl()

```
void Counter::clear_impl ( ) [override], [protected], [virtual]
```

clear [Iterator](#) state.

Each [Iterator](#) should implement the [clear_impl\(\)](#) method to reset its internal state. The [clear_impl\(\)](#) function is guarded by the update lock.

Reimplemented from [_Iterator](#).

7.8.3.2 GET_DATA_2D()

```
Counter::GET_DATA_2D (
    getData ,
    int ,
    array_out )
```

get counts

the counts are copied to a newly allocated allocated memory, an the pointer to this location is returned.

7.8.3.3 next_impl()

```
bool Counter::next_impl (
    std::vector< Tag > & incoming_tags,
    timestamp_t begin_time,
    timestamp_t end_time ) [override], [protected], [virtual]
```

update iterator state

Each [Iterator](#) must implement the [next_impl\(\)](#) method. The [next_impl\(\)](#) function is guarded by the update lock.

The backend delivers each [Tag](#) on each registered channel to this callback function.

Parameters

<i>list</i>	block of events
<i>start_time</i>	earliest event in the block
<i>end_time</i>	start_time of the next block, not including in this block

Returns

if the content of this block was modified

Implements [_Iterator](#).

7.8.3.4 start()

```
void Counter::start ( ) [override], [virtual]
```

start the iterator

The default behavior for iterators is to start automatically on creation.

Reimplemented from [_Iterator](#).

The documentation for this class was generated from the following files:


- Iterators.h
- Counter.cpp

7.9 Countrate Class Reference

count rate on one or more channels

```
#include <Iterators.h>
```

Inheritance diagram for Countrate:



Public Member Functions

- [Countrate](#) ([TimeTagger](#) *tagger, std::vector< channel_t > channels)
constructor of Countrate
- void [start](#) () override
start the iterator
- [GET_DATA_1D](#) (getData, double, array_out,)
get the count rates

Protected Member Functions

- bool `next_impl` (std::vector< [Tag](#) > &incoming_tags, timestamp_t begin_time, timestamp_t end_time) override
update iterator state
- void `clear_impl` () override
clear [Iterator](#) state.

Additional Inherited Members

7.9.1 Detailed Description

count rate on one or more channels

Measures the average count rate on one or more channels. Specifically, it counts incoming clicks and determines the time between the initial click and the latest click. The number of clicks divided by the time corresponds to the average countrate since the initial click.

7.9.2 Constructor & Destructor Documentation

7.9.2.1 Countrate()

```
Countrate::Countrate (
    TimeTagger * tagger,
    std::vector< channel_t > channels )
```

constructor of [Countrate](#)

Parameters

<i>tagger</i>	reference to a TimeTagger
<i>channels</i>	the channels to count on

7.9.3 Member Function Documentation

7.9.3.1 clear_impl()

```
void Countrate::clear_impl ( ) [override], [protected], [virtual]
```

clear [Iterator](#) state.

Each [Iterator](#) should implement the `clear_impl()` method to reset its internal state. The `clear_impl()` function is guarded by the update lock.

Reimplemented from [_Iterator](#).

7.9.3.2 GET_DATA_1D()

```
Countrate::GET_DATA_1D (
    getData ,
    double ,
    array_out )
```

get the count rates

the count rates are copied to a newly allocated allocated memory, an the pointer to this location is returned.

7.9.3.3 next_impl()

```
bool Countrate::next_impl (
    std::vector< Tag > & incoming_tags,
    timestamp_t begin_time,
    timestamp_t end_time ) [override], [protected], [virtual]
```

update iterator state

Each [Iterator](#) must implement the [next_impl\(\)](#) method. The [next_impl\(\)](#) function is guarded by the update lock.

The backend delivers each [Tag](#) on each registered channel to this callback function.

Parameters

<i>list</i>	block of events
<i>start_time</i>	earliest event in the block
<i>end_time</i>	start_time of the next block, not including in this block

Returns

if the content of this block was modified

Implements [_Iterator](#).

7.9.3.4 start()

```
void Countrate::start ( ) [override], [virtual]
```

start the iterator

The default behavior for iterators is to start automatically on creation.

Reimplemented from [_Iterator](#).

The documentation for this class was generated from the following files:


- Iterators.h
- Counter.cpp

7.10 DelayedChannel Class Reference

a simple delayed queue

```
#include <Iterators.h>
```

Inheritance diagram for DelayedChannel:



Public Member Functions

- [DelayedChannel](#) ([TimeTagger](#) *tagger, channel_t input_channel, unsigned long long delay)
constructor of a [DelayedChannel](#)
- void [start](#) () override
start the iterator
- channel_t [getChannel](#) ()
the new virtual channel
- void [setDelay](#) (unsigned long long delay)
set the delay time delay for the cloned tags in the virtual channel $0 \leq t < 2^{63}$

Protected Member Functions

- bool [next_impl](#) (std::vector< [Tag](#) > &incoming_tags, timestamp_t begin_time, timestamp_t end_time) override
update iterator state

Additional Inherited Members

7.10.1 Detailed Description

a simple delayed queue

A simple first-in first-out queue of delayed event timestamps.

7.10.2 Constructor & Destructor Documentation

7.10.2.1 DelayedChannel()

```
DelayedChannel::DelayedChannel (
    TimeTagger * tagger,
    channel_t input_channel,
    unsigned long long delay )
```

constructor of a [DelayedChannel](#)

Parameters

<i>tagger</i>	reference to a TimeTagger
<i>input_channel</i>	channel which is delayed
<i>delay</i>	amount of time to delay, must be positive

7.10.3 Member Function Documentation

7.10.3.1 `getChannel()`

```
channel_t DelayedChannel::getChannel ( )
```

the new virtual channel

This function returns the new allocated virtual channel. It can be used now in any new iterator.

7.10.3.2 `next_impl()`

```
bool DelayedChannel::next_impl (
    std::vector< Tag > & incoming_tags,
    timestamp_t begin_time,
    timestamp_t end_time ) [override], [protected], [virtual]
```

update iterator state

Each [Iterator](#) must implement the `next_impl()` method. The `next_impl()` function is guarded by the update lock.

The backend delivers each [Tag](#) on each registered channel to this callback function.

Parameters

<i>list</i>	block of events
<i>start_time</i>	earliest event in the block
<i>end_time</i>	start_time of the next block, not including in this block

Returns

if the content of this block was modified

Implements [_Iterator](#).

7.10.3.3 setDelay()

```
void DelayedChannel::setDelay (
    unsigned long long delay )
```

set the delay time delay for the cloned tags in the virtual channel $0 \leq t < 2^{63}$

Note: When the delay is the same or greater than the previous value all incoming tags will be visible at virtual channel. By applying a shorter delay time, the tags stored in the local buffer will be flushed and won't be visible in the virtual channel.

7.10.3.4 start()

```
void DelayedChannel::start ( ) [override], [virtual]
```

start the iterator

The default behavior for iterators is to start automatically on creation.

Reimplemented from [_Iterator](#).

The documentation for this class was generated from the following files:


- Iterators.h
- Combiner.cpp

7.11 Dump Class Reference

dump all time tags to a file

```
#include <Iterators.h>
```

Inheritance diagram for Dump:



Public Member Functions

- [Dump](#) ([TimeTagger](#) *tagger, std::string filename, size_t max_tags, const std::vector< channel_t > &channels=std::vector< channel_t >())
constructor of a [Dump](#) thread
- [~Dump](#) ()
tbd

Protected Member Functions

- bool `next_impl` (std::vector< [Tag](#) > &incoming_tags, timestamp_t begin_time, timestamp_t end_time) override
update iterator state
- void `clear_impl` () override
clear [Iterator](#) state.

Additional Inherited Members

7.11.1 Detailed Description

dump all time tags to a file

7.11.2 Constructor & Destructor Documentation

7.11.2.1 ~Dump()

```
Dump::~Dump ( )
```

tbd

7.11.3 Member Function Documentation

7.11.3.1 clear_impl()

```
void Dump::clear_impl ( ) [override], [protected], [virtual]
```

clear [Iterator](#) state.

Each [Iterator](#) should implement the `clear_impl()` method to reset its internal state. The `clear_impl()` function is guarded by the update lock.

Reimplemented from [_Iterator](#).

7.11.3.2 next_impl()

```
bool Dump::next_impl (
    std::vector< Tag > & incoming_tags,
    timestamp_t begin_time,
    timestamp_t end_time ) [override], [protected], [virtual]
```

update iterator state

Each [Iterator](#) must implement the `next_impl()` method. The `next_impl()` function is guarded by the update lock.

The backend delivers each [Tag](#) on each registered channel to this callback function.

Parameters

<i>list</i>	block of events
<i>start_time</i>	earliest event in the block
<i>end_time</i>	start_time of the next block, not including in this block

Returns

if the content of this block was modified

Implements [_iterator](#).

The documentation for this class was generated from the following files:

- Iterators.h
- Iterator.cpp

7.12 Event Struct Reference

Public Attributes

- timestamp_t **time**
- State **state**

The documentation for this struct was generated from the following file:

- Iterators.h

7.13 Flim Class Reference

Fluorescence lifetime imaging.

```
#include <Iterators.h>
```

Public Member Functions

- [Flim](#) ([TimeTagger](#) *tagger, channel_t click_channel, channel_t start_channel, channel_t next_channel, timestamp_t binwidth=1000, int n_bins=1000, int n_pixels=1)
constructor of a FLIM measurement
- void **start** ()
- void **startFor** (timestamp_t capture_duration, bool clear=true)
- void **stop** ()
- void **clear** ()
- timestamp_t **getCaptureDuration** ()
- **GET_DATA_2D** (getData, int, array_out,)
- **GET_DATA_1D** (getIndex, timestamp_t, array_out,)

7.13.1 Detailed Description

Fluorescence lifetime imaging.

Successively acquires n histograms (one for each pixel in the image), where each histogram is determined by the number of bins and the binwidth. Clicks that fall outside the histogram range are ignored.

Fluorescence-lifetime imaging microscopy or FLIM is an imaging technique for producing an image based on the differences in the exponential decay rate of the fluorescence from a fluorescent sample.

Fluorescence lifetimes can be determined in the time domain by using a pulsed source. When a population of fluorophores is excited by an ultrashort or delta pulse of light, the time-resolved fluorescence will decay exponentially.

7.13.2 Constructor & Destructor Documentation

7.13.2.1 Flim()

```
Flim::Flim (
    TimeTagger * tagger,
    channel_t click_channel,
    channel_t start_channel,
    channel_t next_channel,
    timestamp_t binwidth = 1000,
    int n_bins = 1000,
    int n_pixels = 1 )
```

constructor of a FLIM measurement

Parameters

<i>tagger</i>	reference to a TimeTagger
<i>click_channel</i>	channel that increments the count in a bin
<i>start_channel</i>	channel that sets start times relative to which clicks on the click channel are measured
<i>next_channel</i>	channel that increments the pixel
<i>binwidth</i>	width of one histogram bin in ps
<i>n_bins</i>	number of bins in each histogram
<i>n_pixels</i>	number of pixels

The documentation for this class was generated from the following files:


- Iterators.h
- TimeDifferences.cpp

7.14 FrequencyMultiplier Class Reference

The signal of an input channel is scaled up to a higher frequency according to the multiplier passed as a parameter.

```
#include <Iterators.h>
```

Inheritance diagram for FrequencyMultiplier:



Public Member Functions

- [FrequencyMultiplier](#) ([TimeTagger](#) *tagger, channel_t input_channel, int multiplier)
constructor of a [FrequencyMultiplier](#)
- channel_t **getChannel** ()
- int **getMultiplier** ()

Protected Member Functions

- bool [next_impl](#) (std::vector< [Tag](#) > &incoming_tags, timestamp_t begin_time, timestamp_t end_time) override
update iterator state

Additional Inherited Members

7.14.1 Detailed Description

The signal of an input channel is scaled up to a higher frequency according to the multiplier passed as a parameter.

The [FrequencyMultiplier](#) inserts copies the original input events from the input_channel and adds additional events to match the upscaling factor. The algorithm used assumes a constant frequency and calculates out of the last two incoming events linearly the intermediate timestamps to match the upscaled frequency given by the multiplier parameter.

The [FrequencyMultiplier](#) can be used to restore the actual frequency applied to an input_channel which was reduces via the EventDivider to lower the effective data rate. For example a 80 MHz laser sync signal can be scaled down via setEventDivider(..., 80) to 1 MHz (hardware side) and an 80 MHz signal can be restored via [FrequencyMultiplier](#)(..., 80) on the software side with some loss in precision. The [FrequencyMultiplier](#) is an alternative way to reduce the data rate in comparison to the EventFilter, which has a higher precision but can be more difficult to use.

7.14.2 Constructor & Destructor Documentation

7.14.2.1 FrequencyMultiplier()

```
FrequencyMultiplier::FrequencyMultiplier (
    TimeTagger * tagger,
    channel_t input_channel,
    int multiplier )
```

constructor of a [FrequencyMultiplier](#)

Parameters

<i>tagger</i>	reference to a TimeTagger
<i>input_channel</i>	channel on which the upscaling of the frequency is based on
<i>multiplier</i>	frequency upscaling factor

7.14.3 Member Function Documentation

7.14.3.1 next_impl()

```
bool FrequencyMultiplier::next_impl (
    std::vector< Tag > & incoming_tags,
    timestamp_t begin_time,
    timestamp_t end_time ) [override], [protected], [virtual]
```

update iterator state

Each [Iterator](#) must implement the [next_impl\(\)](#) method. The [next_impl\(\)](#) function is guarded by the update lock.

The backend delivers each [Tag](#) on each registered channel to this callback function.

Parameters

<i>list</i>	block of events
<i>start_time</i>	earliest event in the block
<i>end_time</i>	start_time of the next block, not including in this block

Returns

if the content of this block was modified

Implements [_Iterator](#).

The documentation for this class was generated from the following files:


- Iterators.h
- Combiner.cpp

7.15 GatedChannel Class Reference

An input channel is gated by a gate channel.

```
#include <Iterators.h>
```

Inheritance diagram for GatedChannel:



Public Member Functions

- [GatedChannel](#) ([TimeTagger](#) *tagger, channel_t input_channel, channel_t gate_start_channel, channel_t gate_stop_channel)
constructor of a [GatedChannel](#)
- channel_t [getChannel](#) ()
the new virtual channel

Protected Member Functions

- bool [next_impl](#) (std::vector< [Tag](#) > &incoming_tags, timestamp_t begin_time, timestamp_t end_time) override
update iterator state

Additional Inherited Members

7.15.1 Detailed Description

An input channel is gated by a gate channel.

Note: The gate is edge sensitive and not level sensitive. That means that the gate will transfer data only when an appropriate level change is detected on the gate_start_channel.

7.15.2 Constructor & Destructor Documentation

7.15.2.1 GatedChannel()

```

GatedChannel::GatedChannel (
    TimeTagger * tagger,
    channel_t input_channel,
    channel_t gate_start_channel,
    channel_t gate_stop_channel )
  
```

constructor of a [GatedChannel](#)

Parameters

<i>tagger</i>	reference to a TimeTagger
<i>input_channel</i>	channel which is gated
<i>gate_start_channel</i>	channel on which a signal detected will start the transmission of the input_channel through the gate
<i>gate_stop_channel</i>	channel on which a signal detected will stop the transmission of the input_channel through the gate

7.15.3 Member Function Documentation

7.15.3.1 `getChannel()`

```
channel_t GatedChannel::getChannel ( )
```

the new virtual channel

This function returns the new allocated virtual channel. It can be used now in any new iterator.

7.15.3.2 `next_impl()`

```
bool GatedChannel::next_impl (
    std::vector< Tag > & incoming_tags,
    timestamp_t begin_time,
    timestamp_t end_time ) [override], [protected], [virtual]
```

update iterator state

Each [Iterator](#) must implement the `next_impl()` method. The `next_impl()` function is guarded by the update lock.

The backend delivers each [Tag](#) on each registered channel to this callback function.

Parameters

<i>list</i>	block of events
<i>start_time</i>	earliest event in the block
<i>end_time</i>	start_time of the next block, not including in this block

Returns

if the content of this block was modified

Implements [_Iterator](#).

The documentation for this class was generated from the following files:

- `Iterators.h`
- `Combiner.cpp`

7.16 Histogram Class Reference

Accumulate time differences into a histogram.

```
#include <Iterators.h>
```

Public Member Functions

- [Histogram](#) ([TimeTagger](#) *tagger, channel_t click_channel, channel_t start_channel=CHANNEL_UNUSED, timestamp_t binwidth=1000, int n_bins=1000)
constructor of a [Histogram](#) measurement
- void **start** ()
- void **startFor** (timestamp_t capture_duration, bool clear=true)
- void **stop** ()
- void **clear** ()
- timestamp_t **getCaptureDuration** ()
- bool **isRunning** ()
- **GET_DATA_1D** (getData, int, array_out,)
- **GET_DATA_1D** (getIndex, timestamp_t, array_out,)

7.16.1 Detailed Description

Accumulate time differences into a histogram.

This is a simple multiple start, multiple stop measurement. This is a special case of the more general '[TimeDifferences](#)' measurement. Specifically, the thread waits for clicks on a first channel, the 'start channel', then measures the time difference between the last start click and all subsequent clicks on a second channel, the 'click channel', and stores them in a histogram. The histogram range and resolution is specified by the number of bins and the binwidth. Clicks that fall outside the histogram range are ignored. Data accumulation is performed independently for all start clicks. This type of measurement is frequently referred to as 'multiple start, multiple stop' measurement and corresponds to a full auto- or cross-correlation measurement.

7.16.2 Constructor & Destructor Documentation

7.16.2.1 Histogram()

```
Histogram::Histogram (
    TimeTagger * tagger,
    channel_t click_channel,
    channel_t start_channel = CHANNEL_UNUSED,
    timestamp_t binwidth = 1000,
    int n_bins = 1000 )
```

constructor of a [Histogram](#) measurement

Parameters

<i>tagger</i>	reference to a TimeTagger
<i>click_channel</i>	channel that increments the count in a bin
<i>start_channel</i>	channel that sets start times relative to which clicks on the click channel are measured
<i>binwidth</i>	width of one histogram bin in ps
<i>n_bins</i>	number of bins in the histogram

The documentation for this class was generated from the following files:


- Iterators.h
- TimeDifferences.cpp

7.17 HistogramLogBins Class Reference

Accumulate time differences into a histogram with logarithmic increasing bin sizes.

```
#include <Iterators.h>
```

Inheritance diagram for HistogramLogBins:



Public Member Functions

- [HistogramLogBins](#) ([TimeTagger](#) *tagger, channel_t click_channel, channel_t start_channel, double exp_start, double exp_stop, int n_bins)
constructor of a [HistogramLogBins](#) measurement
- **GET_DATA_1D** (getData, unsigned long long, array_out,)
- **GET_DATA_1D** (getDataNormalized, double, array_out,)
- **GET_DATA_1D** (getBinEdges, timestamp_t, array_out,)

Protected Member Functions

- bool [next_impl](#) (std::vector< [Tag](#) > &incoming_tags, timestamp_t begin_time, timestamp_t end_time) override
update iterator state
- void [clear_impl](#) () override
clear [Iterator](#) state.

Additional Inherited Members

7.17.1 Detailed Description

Accumulate time differences into a histogram with logarithmic increasing bin sizes.

This is a multiple start, multiple stop measurement, and works the very same way as the histogram measurement but with logarithmic increasing bin widths. After initializing the measurement (or after an overflow) no data is accumulated in the histogram until the full histogram duration has passed to ensure a balanced count accumulation over the full histogram.

7.17.2 Constructor & Destructor Documentation

7.17.2.1 HistogramLogBins()

```
HistogramLogBins::HistogramLogBins (
    TimeTagger * tagger,
    channel_t click_channel,
    channel_t start_channel,
    double exp_start,
    double exp_stop,
    int n_bins )
```

constructor of a [HistogramLogBins](#) measurement

Parameters

<i>tagger</i>	reference to a TimeTagger
<i>click_channel</i>	channel that increments the count in a bin
<i>start_channel</i>	channel that sets start times relative to which clicks on the click channel are measured
<i>exp_start</i>	exponent for the lowest time differences in the histogram: $10^{\text{exp_start}}$ s, lowest exp_start: -12 => 1ps
<i>exp_stop</i>	exponent for the highest time differences in the histogram: $10^{\text{exp_stop}}$ s
<i>n_bins</i>	total number of bins in the histogram

7.17.3 Member Function Documentation

7.17.3.1 clear_impl()

```
void HistogramLogBins::clear_impl ( ) [override], [protected], [virtual]
```

clear [Iterator](#) state.

Each [Iterator](#) should implement the [clear_impl\(\)](#) method to reset its internal state. The [clear_impl\(\)](#) function is guarded by the update lock.

Reimplemented from [_Iterator](#).

7.17.3.2 next_impl()

```
bool HistogramLogBins::next_impl (
    std::vector< Tag > & incoming_tags,
    timestamp_t begin_time,
    timestamp_t end_time ) [override], [protected], [virtual]
```

update iterator state

Each [Iterator](#) must implement the [next_impl\(\)](#) method. The [next_impl\(\)](#) function is guarded by the update lock.

The backend delivers each [Tag](#) on each registered channel to this callback function.

Parameters

<i>list</i>	block of events
<i>start_time</i>	earliest event in the block
<i>end_time</i>	start_time of the next block, not including in this block

Returns

if the content of this block was modified

Implements [_Iterator](#).

The documentation for this class was generated from the following files:


- Iterators.h
- TimeDifferences.cpp

7.18 Iterator Class Reference

a simple event queue

```
#include <Iterators.h>
```

Inheritance diagram for Iterator:



Public Member Functions

- [Iterator](#) ([TimeTagger](#) *tagger, channel_t channel)
standard constructor
- timestamp_t [next](#) ()
get next timestamp
- int [size](#) ()
get queue size

Protected Member Functions

- bool [next_impl](#) (std::vector< [Tag](#) > &incoming_tags, timestamp_t begin_time, timestamp_t end_time) override
update iterator state
- void [clear_impl](#) () override
clear [Iterator](#) state.

Additional Inherited Members

7.18.1 Detailed Description

a simple event queue

A simple [Iterator](#), just keeping a first-in first-out queue of event timestamps.

7.18.2 Constructor & Destructor Documentation

7.18.2.1 Iterator()

```
Iterator::Iterator (
    TimeTagger * tagger,
    channel_t channel )
```

standard constructor

Parameters

<i>tagger</i>	the backend
<i>channel</i>	the channel to get events from

7.18.3 Member Function Documentation

7.18.3.1 clear_impl()

```
void Iterator::clear_impl ( ) [override], [protected], [virtual]
```

clear [Iterator](#) state.

Each [Iterator](#) should implement the [clear_impl\(\)](#) method to reset its internal state. The [clear_impl\(\)](#) function is guarded by the update lock.

Reimplemented from [_Iterator](#).

7.18.3.2 next()

```
timestamp_t Iterator::next ( )
```

get next timestamp

get the next timestamp from the queue.

7.18.3.3 next_impl()

```
bool Iterator::next_impl (
    std::vector< Tag > & incoming_tags,
    timestamp_t begin_time,
    timestamp_t end_time ) [override], [protected], [virtual]
```

update iterator state

Each [Iterator](#) must implement the [next_impl\(\)](#) method. The [next_impl\(\)](#) function is guarded by the update lock.

The backend delivers each [Tag](#) on each registered channel to this callback function.

Parameters

<i>list</i>	block of events
<i>start_time</i>	earliest event in the block
<i>end_time</i>	start_time of the next block, not including in this block

Returns

if the content of this block was modified

Implements [_Iterator](#).

7.18.3.4 size()

```
int Iterator::size ( )
```

get queue size

The documentation for this class was generated from the following files:

- Iterators.h
- Iterator.cpp

7.19 LED Union Reference

Public Attributes

•

```

struct {
    uint32_t statusR: 1
    uint32_t statusB: 1
    uint32_t statusG: 1
    uint32_t powerR: 1
    uint32_t powerB: 1
    uint32_t powerG: 1
    uint32_t clockR: 1
    uint32_t clockB: 1
    uint32_t clockG: 1
};

```

- uint32_t **ledStatus**

The documentation for this union was generated from the following file:

- TimeTaggerTest.cpp

7.20 PRBS Class Reference

Public Member Functions


- uint32_t **get** ()

The documentation for this class was generated from the following file:

- PRBS.h

7.21 Scope Class Reference

Inheritance diagram for Scope:



Public Member Functions

- [Scope](#) ([TimeTagger](#) *tagger, std::vector< channel_t > event_channels, channel_t trigger_channel, timestamp_t window_size=1000000000, int n_traces=1, int n_max_events=1000)
constructor of a [Scope](#) measurement
- bool **ready** ()
- int **triggered** ()
- void [start](#) () override
start the iterator
- std::vector< std::vector< [Event](#) > > **getData** ()
- timestamp_t **getWindowSize** ()

Protected Member Functions

- bool `next_impl` (std::vector< [Tag](#) > &incoming_tags, timestamp_t begin_time, timestamp_t end_time) override
update iterator state
- void `clear_impl` () override
clear [Iterator](#) state.

Additional Inherited Members

7.21.1 Constructor & Destructor Documentation

7.21.1.1 Scope()

```
Scope::Scope (
    TimeTagger * tagger,
    std::vector< channel_t > event_channels,
    channel_t trigger_channel,
    timestamp_t window_size = 1000000000,
    int n_traces = 1,
    int n_max_events = 1000 )
```

constructor of a [Scope](#) measurement

Parameters

<i>tagger</i>	reference to a TimeTagger
<i>event_channels</i>	channels which are captured
<i>trigger_channel</i>	channel that starts a new trace
<i>window_size</i>	window time of each trace
<i>n_traces</i>	amount of traces ($n_traces < 1$, automatic retrigger)
<i>n_max_events</i>	maximum number of tags in each trace

7.21.2 Member Function Documentation

7.21.2.1 clear_impl()

```
void Scope::clear_impl ( ) [override], [protected], [virtual]
```

clear [Iterator](#) state.

Each [Iterator](#) should implement the `clear_impl()` method to reset its internal state. The `clear_impl()` function is guarded by the update lock.

Reimplemented from [_Iterator](#).

7.21.2.2 next_impl()

```
bool Scope::next_impl (
    std::vector< Tag > & incoming_tags,
    timestamp_t begin_time,
    timestamp_t end_time ) [override], [protected], [virtual]
```

update iterator state

Each [Iterator](#) must implement the [next_impl\(\)](#) method. The [next_impl\(\)](#) function is guarded by the update lock.

The backend delivers each [Tag](#) on each registered channel to this callback function.

Parameters

<i>list</i>	block of events
<i>start_time</i>	earliest event in the block
<i>end_time</i>	start_time of the next block, not including in this block

Returns

if the content of this block was modified

Implements [_Iterator](#).

7.21.2.3 start()

```
void Scope::start ( ) [override], [virtual]
```

start the iterator

The default behavior for iterators is to start automatically on creation.

Reimplemented from [_Iterator](#).

The documentation for this class was generated from the following files:


- Iterators.h
- Scope.cpp

7.22 StartStop Class Reference

simple start-stop measurement

```
#include <Iterators.h>
```

Inheritance diagram for StartStop:



Public Member Functions

- [StartStop](#) ([TimeTagger](#) *tagger, channel_t click_channel, channel_t start_channel=CHANNEL_UNUSED, timestamp_t binwidth=1000)
constructor of [StartStop](#)
- void [start](#) () override
start the iterator
- **GET_DATA_2D** (getData, timestamp_t, array_out,)

Protected Member Functions

- bool [next_impl](#) (std::vector< [Tag](#) > &incoming_tags, timestamp_t begin_time, timestamp_t end_time) override
update iterator state
- void [clear_impl](#) () override
clear [Iterator](#) state.

Additional Inherited Members

7.22.1 Detailed Description

simple start-stop measurement

This class performs a start-stop measurement between two channels and stores the time differences in a histogram. The histogram resolution is specified beforehand (binwidth) but the histogram range is unlimited. It is adapted to the largest time difference that was detected. Thus all pairs of subsequent clicks are registered.

Be aware, on long-running measurements this may considerably slow down system performance and even crash the system entirely when attached to an unsuitable signal source.

7.22.2 Constructor & Destructor Documentation

7.22.2.1 StartStop()

```
StartStop::StartStop (
    TimeTagger * tagger,
    channel_t click_channel,
    channel_t start_channel = CHANNEL_UNUSED,
    timestamp_t binwidth = 1000 )
```

constructor of [StartStop](#)

Parameters

<i>tagger</i>	reference to a TimeTagger
<i>click_channel</i>	channel for stop clicks
<i>start_channel</i>	channel for start clicks
<i>binwidth</i>	width of one histogram bin in ps

7.22.3 Member Function Documentation

7.22.3.1 clear_impl()

```
void StartStop::clear_impl ( ) [override], [protected], [virtual]
```

clear [Iterator](#) state.

Each [Iterator](#) should implement the [clear_impl\(\)](#) method to reset its internal state. The [clear_impl\(\)](#) function is guarded by the update lock.

Reimplemented from [_Iterator](#).

7.22.3.2 next_impl()

```
bool StartStop::next_impl (
    std::vector< Tag > & incoming_tags,
    timestamp_t begin_time,
    timestamp_t end_time ) [override], [protected], [virtual]
```

update iterator state

Each [Iterator](#) must implement the [next_impl\(\)](#) method. The [next_impl\(\)](#) function is guarded by the update lock.

The backend delivers each [Tag](#) on each registered channel to this callback function.

Parameters

<i>list</i>	block of events
<i>start_time</i>	earliest event in the block
<i>end_time</i>	start_time of the next block, not including in this block

Returns

if the content of this block was modified

Implements [_Iterator](#).

7.22.3.3 start()

```
void StartStop::start ( ) [override], [virtual]
```

start the iterator

The default behavior for iterators is to start automatically on creation.

Reimplemented from [_Iterator](#).

The documentation for this class was generated from the following files:

- Iterators.h
- StartStop.cpp

7.23 SynchronizedMeasurements Class Reference

start, stop and clear several measurements synchronized

```
#include <Iterators.h>
```

Public Member Functions

- [SynchronizedMeasurements](#) ([TimeTagger](#) *tagger)
construct a SynchronizedMeasurements object - if you have an better idea how to call it please let me know
- void [registerMeasurement](#) ([_Iterator](#) *measurement)
register a measurement (iterator) to the SynchronizedMeasurements-group.
- void [clear](#) ()
clear all registered measurements synchronously
- void [start](#) ()
start all registered measurements synchronously
- void [stop](#) ()
stop all registered measurements synchronously

7.23.1 Detailed Description

start, stop and clear several measurements synchronized

For the case that several measurements should be started, stopped or cleared at the very same time, a SynchronizedMeasurements object can be create to which all the measurements (also called iterators) can be registered with `.registerMeasurement(measurement)`. Calling `.stop()`, `.start()` or `.clear()` on the SynchronizedMeasurements object will call the respective method on each of the registered measurements at the very same time. That means that all measurements taking part will have processed the very same time tags.

7.23.2 Constructor & Destructor Documentation

7.23.2.1 SynchronizedMeasurements()

```
SynchronizedMeasurements::SynchronizedMeasurements (
    TimeTagger * tagger ) [inline]
```

construct a SynchronizedMeasurements object - if you have an better idea how to call it please let me know

Parameters

<i>tagger</i>	reference to a TimeTagger
---------------	---

7.23.3 Member Function Documentation

7.23.3.1 clear()

```
void SynchronizedMeasurments::clear ( ) [inline]
```

clear all registered measurements synchronously

7.23.3.2 registerMeasurement()

```
void SynchronizedMeasurments::registerMeasurement (
    _Iterator * measurement ) [inline]
```

register a measurement (iterator) to the SynchronizedMeasurements-group.

All available methods called on the SynchronizedMeasurements will happen at the very same time for all the registered measurements.

7.23.3.3 start()

```
void SynchronizedMeasurments::start ( ) [inline]
```

start all registered measurements synchronously

7.23.3.4 stop()

```
void SynchronizedMeasurments::stop ( ) [inline]
```

stop all registered measurements synchronously

The documentation for this class was generated from the following file:

- Iterators.h

7.24 Tag Struct Reference

a single event on a channel

```
#include <TimeTagger.h>
```

Public Attributes

- bool [overflow](#)
when set, there was an overflow on the communication channel.
- channel_t [channel](#)
the channel number
- timestamp_t [time](#)
the timestamp on the event, in picoseconds

7.24.1 Detailed Description

a single event on a channel

Channel events are passed from the backend to registered iterators by the `_Iterator::next()` callback function.

A [Tag](#) describes a single event on a channel.

7.24.2 Member Data Documentation

7.24.2.1 channel

```
channel_t Tag::channel
```

the channel number

7.24.2.2 overflow

```
bool Tag::overflow
```

when set, there was an overflow on the communication channel.

7.24.2.3 time

```
timestamp_t Tag::time
```

the timestamp on the event, in picoseconds

The documentation for this struct was generated from the following file:


- TimeTagger.h

7.25 TimeDifferences Class Reference

Accumulates the time differences between clicks on two channels in one or more histograms.

```
#include <Iterators.h>
```

Inheritance diagram for TimeDifferences:



Public Member Functions

- **TimeDifferences** ([TimeTagger](#) *tagger, channel_t click_channel, channel_t start_channel=CHANNEL_UNUSED, channel_t next_channel=CHANNEL_UNUSED, channel_t sync_channel=CHANNEL_UNUSED, timestamp_t binwidth=1000, int n_bins=1000, int n_histograms=1)
constructor of a [TimeDifferences](#) measurement
- void [start](#) () override
start the iterator
- [GET_DATA_2D](#) (getData, int, array_out,)
get result data
- **GET_DATA_1D** (getIndex, timestamp_t, array_out,)
- void [setMaxCounts](#) (int max_counts)
set the number of sync/next clicks after which acquisition shall stop
- int [getCounts](#) ()
return the number of sync/next clicks
- bool [ready](#) ()
return 'true' if the maximum number of start clicks has been reached

Protected Member Functions

- bool [next_impl](#) (std::vector< [Tag](#) > &incoming_tags, timestamp_t begin_time, timestamp_t end_time) override
update iterator state
- void [clear_impl](#) () override
clear [Iterator](#) state.

Additional Inherited Members

7.25.1 Detailed Description

Accumulates the time differences between clicks on two channels in one or more histograms.

Specifically, the thread waits for clicks on a first channel, the 'start channel', then measures the time difference between a start click and all subsequent clicks on a second channel, the 'click channel', and stores them in a histogram. The histogram range and resolution is specified by the number of bins and the binwidth. Clicks that fall outside the histogram range are ignored. Data accumulation is performed independently for all start clicks. This type of measurement is frequently referred to as 'multiple start, multiple stop' measurement and corresponds to a full auto- or cross-correlation measurement. The data obtained from subsequent start clicks can be accumulated into the same histogram (one-dimensional measurement) or into different histograms (two-dimensional measurement). In this way you can perform more complex two-dimensional time-difference measurements. Specifically, after each click on the next channel, the histogram index is incremented by one and reset to zero after the last histogram has been served. You can also provide an external synchronization trigger that resets the histogram index to zero.

7.25.2 Constructor & Destructor Documentation

7.25.2.1 TimeDifferences()

```
TimeDifferences::TimeDifferences (
    TimeTagger * tagger,
    channel_t click_channel,
    channel_t start_channel = CHANNEL_UNUSED,
    channel_t next_channel = CHANNEL_UNUSED,
    channel_t sync_channel = CHANNEL_UNUSED,
    timestamp_t binwidth = 1000,
    int n_bins = 1000,
    int n_histograms = 1 )
```

constructor of a [TimeDifferences](#) measurement

Parameters

<i>tagger</i>	reference to a TimeTagger
<i>click_channel</i>	channel that increments the count in a bin
<i>start_channel</i>	channel that sets start times relative to which clicks on the click channel are measured
<i>next_channel</i>	channel that increments the histogram index
<i>sync_channel</i>	channel that resets the histogram index to zero
<i>binwidth</i>	width of one histogram bin in ps
<i>n_bins</i>	number of bins in each histogram
<i>n_histograms</i>	number of histograms

7.25.3 Member Function Documentation

7.25.3.1 clear_impl()

```
void TimeDifferences::clear_impl ( ) [override], [protected], [virtual]
```

clear [Iterator](#) state.

Each [Iterator](#) should implement the [clear_impl\(\)](#) method to reset its internal state. The [clear_impl\(\)](#) function is guarded by the update lock.

Reimplemented from [_Iterator](#).

7.25.3.2 GET_DATA_2D()

```
TimeDifferences::GET_DATA_2D (
    getData ,
    int ,
    array_out )
```

get result data

Parameters

<i>ARGOUTVIEWM_ARRAY2</i>	pointer receiving pointer to data
<i>DIM1</i>	pointer receiving first dimension
<i>DIM2</i>	pointer receiving second dimension

7.25.3.3 getCounts()

```
int TimeDifferences::getCounts ( )
```

return the number of sync/next clicks

7.25.3.4 next_impl()

```
bool TimeDifferences::next_impl (
    std::vector< Tag > & incoming_tags,
    timestamp_t begin_time,
    timestamp_t end_time ) [override], [protected], [virtual]
```

update iterator state

Each [Iterator](#) must implement the [next_impl\(\)](#) method. The [next_impl\(\)](#) function is guarded by the update lock.

The backend delivers each [Tag](#) on each registered channel to this callback function.

Parameters

<i>list</i>	block of events
<i>start_time</i>	earliest event in the block
<i>end_time</i>	start_time of the next block, not including in this block

Returns

if the content of this block was modified

Implements [_Iterator](#).

7.25.3.5 ready()

```
bool TimeDifferences::ready ( )
```

return 'true' if the maximum number of start clicks has been reached

7.25.3.6 setMaxCounts()

```
void TimeDifferences::setMaxCounts (
    int max_counts )
```

set the number of sync/next clicks after which acquisition shall stop

Parameters

<i>counts</i>	maximum number of sync/next clicks
---------------	------------------------------------

7.25.3.7 start()

```
void TimeDifferences::start ( ) [override], [virtual]
```

start the iterator

The default behavior for iterators is to start automatically on creation.

Reimplemented from [_Iterator](#).

The documentation for this class was generated from the following files:

- Iterators.h
- TimeDifferences.cpp

7.26 TimeTagger Class Reference

backend for the [TimeTagger](#).

```
#include <TimeTagger.h>
```

Public Types

- typedef std::function< void([_Iterator](#) *)> **IteratorCallback**
- typedef std::map< [_Iterator](#) *, IteratorCallback > **IteratorCallbackMap**

Public Member Functions

- virtual void [reset](#) ()=0
reset the [TimeTagger](#) object to default settings and detach all iterators
- virtual void [setTestSignalDivider](#) (int divider)=0
set the divider for the frequency of the test signal
- virtual void [setTriggerLevel](#) (channel_t channel, double voltage)=0
set the trigger voltage threshold of a channel
- virtual double [getTriggerLevel](#) (channel_t channel)=0
get the trigger voltage threshold of a channel
- virtual void [setInputDelay](#) (channel_t channel, timestamp_t delay)=0
set time delay on a channel
- virtual timestamp_t [getInputDelay](#) (channel_t channel)=0
get time delay of a channel
- virtual void [setConditionalFilter](#) (std::vector< channel_t > trigger, std::vector< channel_t > filtered)=0
configures the conditional filter
- virtual std::vector< channel_t > [getConditionalFilterTrigger](#) ()=0
fetches the configuration of the conditional filter
- virtual std::vector< channel_t > [getConditionalFilterFiltered](#) ()=0
fetches the configuration of the conditional filter
- virtual void [setFilter](#) (bool state)=0
enables or disables the filter on the FPGA board.
- virtual bool [getFilter](#) ()=0
returns the filter state on the FPGA board
- virtual void [setNormalization](#) (bool state)=0
enables or disables the normalization of the distribution.
- virtual bool [getNormalization](#) ()=0
returns the the normalization of the distribution.
- virtual void [setHardwareBufferSize](#) (int size)=0
sets the maximum USB buffer size
- virtual int [getHardwareBufferSize](#) ()=0
queries the size of the USB queue
- virtual timestamp_t [setDeadtime](#) (channel_t channel, timestamp_t deadtime)=0
set the deadtime between two edges on the same channel.
- virtual timestamp_t [getDeadtime](#) (channel_t channel)=0
get the deadtime between two edges on the same channel.
- virtual void [setEventDivider](#) (channel_t channel, unsigned int divider)=0
Divides the amount of transmitted edge per channel.

- virtual unsigned int [getEventDivider](#) (channel_t channel)=0
Returns the factor of the dividing filter.
- virtual void [registerChannel](#) (channel_t channel)=0
register a FPGA channel.
- virtual void [unregisterChannel](#) (channel_t channel)=0
release a previously registered channel.
- virtual void [setTestSignal](#) (channel_t channel, bool enabled)=0
enable the calibration on a channel.
- virtual void [setTestSignal](#) (std::vector< channel_t > channel, bool enabled)=0
- virtual bool [getTestSignal](#) (channel_t channel)=0
fetch the status of the test signal generator
- virtual void [autoCalibration](#) (bool verbose=true)=0
runs a calibrations based on the on-chip uncorrelated signal generator.
- virtual std::string [getSerial](#) ()=0
identifies the hardware by serial number
- virtual std::string [getModel](#) ()=0
identifies the hardware by Time Tagger Model
- virtual int [getChannelNumberScheme](#) ()=0
Fetch the configured numbering scheme for this TimeTagger object.
- virtual std::vector< double > [getDACRange](#) ()=0
returns the minimum and the maximum voltage of the DACs as a trigger reference
- [GET_DATA_2D](#) (getDistributionCount, long long, array_out, virtual)=0
get internal calibration data
- [GET_DATA_2D](#) (getDistributionPSecs, timestamp_t, array_out, virtual)=0
get internal calibration data
- virtual channel_t [getChannels](#) ()=0
fetch the amount of channels
- virtual std::vector< channel_t > [getChannelList](#) (int type=TT_CHANNEL_RISING_AND_FALLING_EDGE←GES)=0
- virtual channel_t [getInvertedChannel](#) (channel_t channel)=0
get the falling channel id for a raising channel and vice versa
- virtual bool [isUnusedChannel](#) (channel_t channel)=0
compares the provided channel with CHANNEL_UNUSED
- virtual timestamp_t [getPsPerClock](#) ()=0
fetch the duration of each clock cycle in picoseconds
- virtual long long [getOverflows](#) ()=0
get overflow count
- virtual void [clearOverflows](#) ()=0
clear overflow counter
- virtual long long [getOverflowsAndClear](#) ()=0
get and clear overflow counter
- virtual void [sync](#) ()=0
Sync the timetagger pipeline, so that all started iterators and their enabled channels are ready.
- virtual std::string [getPcbVersion](#) ()=0
Return the hardware version of the PCB board. Version 0 is everything before mid 2018 and with the channel configuration ZERO. version >= 1 is channel configuration ONE.
- virtual std::string [getSensorData](#) ()=0
Show the status of the sensor data from the FPGA and peripherals on the console.
- virtual void [setLED](#) (uint32_t bitmask)=0

Enforce a state to the LEDs 0: `led_status[R]` 16: `led_status[R]` - mux 1: `led_status[B]` 17: `led_status[B]` - mux 2: `led_status[G]` 18: `led_status[G]` - mux 3: `led_power[R]` 19: `led_power[R]` - mux 4: `led_power[B]` 20: `led_power[B]` - mux 5: `led_power[G]` 21: `led_power[G]` - mux 6: `led_clock[R]` 22: `led_clock[R]` - mux 7: `led_clock[B]` 23: `led_clock[B]` - mux 8: `led_clock[G]` 24: `led_clock[G]` - mux.

- virtual void `runSynchronized` (const `IteratorCallbackMap` &callbacks, bool block=true)=0

Run synchronized callbacks for a list of iterators.

Protected Member Functions

- `TimeTagger` ()
abstract interface class
- virtual `~TimeTagger` ()
- `TimeTagger` (const `TimeTagger` &)=delete
- `TimeTagger` & `operator=` (const `TimeTagger` &)=delete

Friends

- class `_Iterator`

7.26.1 Detailed Description

backend for the `TimeTagger`.

The `TimeTagger` class connects to the hardware, and handles the communication over the usb. There may be only one instance of the backend per physical device.

7.26.2 Constructor & Destructor Documentation

7.26.2.1 `TimeTagger()`

```
TimeTagger::TimeTagger ( ) [inline], [protected]
```

abstract interface class

7.26.2.2 `~TimeTagger()`

```
virtual TimeTagger::~~TimeTagger ( ) [inline], [protected], [virtual]
```

destructor

7.26.3 Member Function Documentation

7.26.3.1 `autoCalibration()`

```
virtual void TimeTagger::autoCalibration (
    bool verbose = true ) [pure virtual]
```

runs a calibrations based on the on-chip uncorrelated signal generator.

Parameters

<i>verbose</i>	Verbose output on stdout
----------------	--------------------------

7.26.3.2 clearOverflows()

```
virtual void TimeTagger::clearOverflows ( ) [pure virtual]
```

clear overflow counter

Sets the overflow counter to zero

7.26.3.3 GET_DATA_2D() [1/2]

```
TimeTagger::GET_DATA_2D (
    getDistributionCount ,
    long long,
    array_out ,
    virtual ) [pure virtual]
```

get internal calibration data

7.26.3.4 GET_DATA_2D() [2/2]

```
TimeTagger::GET_DATA_2D (
    getDistributionPSecs ,
    timestamp_t ,
    array_out ,
    virtual ) [pure virtual]
```

get internal calibration data

7.26.3.5 getChannelNumberScheme()

```
virtual int TimeTagger::getChannelNumberScheme ( ) [pure virtual]
```

Fetch the configured numbering scheme for this [TimeTagger](#) object.

Please see setTimeTaggerChannelNumberScheme() for details.

7.26.3.6 getChannels()

```
virtual channel_t TimeTagger::getChannels ( ) [pure virtual]
```

fetch the amount of channels

Deprecated Use getChannelList instead.

7.26.3.7 getConditionalFilterFiltered()

```
virtual std::vector<channel_t> TimeTagger::getConditionalFilterFiltered ( ) [pure virtual]
```

fetches the configuration of the conditional filter

see setConditionalFilter

7.26.3.8 getConditionalFilterTrigger()

```
virtual std::vector<channel_t> TimeTagger::getConditionalFilterTrigger ( ) [pure virtual]
```

fetches the configuration of the conditional filter

see setConditionalFilter

7.26.3.9 getDACRange()

```
virtual std::vector<double> TimeTagger::getDACRange ( ) [pure virtual]
```

returns the minumum and the maximum voltage of the DACs as a trigger reference

7.26.3.10 getDeadtime()

```
virtual timestamp_t TimeTagger::getDeadtime (
    channel_t channel ) [pure virtual]
```

get the deadtime between two edges on the same channel.

This function gets the user configureable deadtime.

Parameters

<i>channel</i>	channel to be queried
----------------	-----------------------

Returns

the real configured deadtime

7.26.3.11 getEventDivider()

```
virtual unsigned int TimeTagger::getEventDivider (
    channel_t channel ) [pure virtual]
```

Returns the factor of the dividing filter.

See for further details.

Parameters

<i>channel</i>	channel to be queried
----------------	-----------------------

Returns

the configured divider

7.26.3.12 getFilter()

```
virtual bool TimeTagger::getFilter ( ) [pure virtual]
```

returns the filter state on the FPGA board

Deprecated use getConditionalFilter*

The laserfilter disables transmission of nontriggered tags on channel 7. This is a deprecated specialization of getConditionalFilter*.

7.26.3.13 getHardwareBufferSize()

```
virtual int TimeTagger::getHardwareBufferSize ( ) [pure virtual]
```

queries the size of the USB queue

See for more information.

Returns

the actual size of the USB queue in events

7.26.3.14 getInputDelay()

```
virtual timestamp_t TimeTagger::getInputDelay (
    channel_t channel ) [pure virtual]
```

get time delay of a channel

see setInputDelay

Parameters

<i>channel</i>	the channel
----------------	-------------

7.26.3.15 getInvertedChannel()

```
virtual channel_t TimeTagger::getInvertedChannel (
    channel_t channel ) [pure virtual]
```

get the falling channel id for a raising channel and vice versa

7.26.3.16 getModel()

```
virtual std::string TimeTagger::getModel ( ) [pure virtual]
```

identifies the hardware by Time Tagger Model

7.26.3.17 getNormalization()

```
virtual bool TimeTagger::getNormalization ( ) [pure virtual]
```

returns the the normalization of the distribution.

Refer the Manual for a description of this function.

7.26.3.18 getOverflows()

```
virtual long long TimeTagger::getOverflows ( ) [pure virtual]
```

get overflow count

Get the number of communication overflows occurred

7.26.3.19 getOverflowsAndClear()

```
virtual long long TimeTagger::getOverflowsAndClear ( ) [pure virtual]
```

get and clear overflow counter

Get the number of communication overflows occurred and sets them to zero

7.26.3.20 getPcbVersion()

```
virtual std::string TimeTagger::getPcbVersion ( ) [pure virtual]
```

Return the hardware version of the PCB board. Version 0 is everything before mid 2018 and with the channel configuration ZERO. version ≥ 1 is channel configuration ONE.

7.26.3.21 getPsPerClock()

```
virtual timestamp_t TimeTagger::getPsPerClock ( ) [pure virtual]
```

fetch the duration of each clock cycle in picoseconds

7.26.3.22 getSensorData()

```
virtual std::string TimeTagger::getSensorData ( ) [pure virtual]
```

Show the status of the sensor data from the FPGA and peripherals on the console.

7.26.3.23 getSerial()

```
virtual std::string TimeTagger::getSerial ( ) [pure virtual]
```

identifies the hardware by serial number

7.26.3.24 getTestSignal()

```
virtual bool TimeTagger::getTestSignal (
    channel_t channel ) [pure virtual]
```

fetch the status of the test signal generator

Parameters

<i>channel</i>	the channel
----------------	-------------

7.26.3.25 getTriggerLevel()

```
virtual double TimeTagger::getTriggerLevel (
    channel_t channel ) [pure virtual]
```

get the trigger voltage threshold of a channel

Parameters

<i>channel</i>	the channel
----------------	-------------

7.26.3.26 isUnusedChannel()

```
virtual bool TimeTagger::isUnusedChannel (
    channel_t channel ) [pure virtual]
```

compares the provided channel with CHANNEL_UNUSED

But also keeps care about the channel number scheme and selects either CHANNEL_UNUSED or CHANNEL_UNUSED_OLD

7.26.3.27 registerChannel()

```
virtual void TimeTagger::registerChannel (
    channel_t channel ) [pure virtual]
```

register a FPGA channel.

Only events on previously registered channels will be transfered over the communication channel.

Parameters

<i>channel</i>	the channel
----------------	-------------

7.26.3.28 reset()

```
virtual void TimeTagger::reset ( ) [pure virtual]
```

reset the [TimeTagger](#) object to default settings and detach all iterators

7.26.3.29 runSynchronized()

```
virtual void TimeTagger::runSynchronized (
    const IteratorCallbackMap & callbacks,
    bool block = true ) [pure virtual]
```

Run synchronized callbacks for a list of iterators.

This method has a list of callbacks for a list of iterators. Those callbacks are called for a synchronized data set, but in parallel. They are called from an internal worker thread. As the data set is synchronized, this creates a bottleneck for one worker thread, so only fast and non-blocking callbacks are allowed.

Parameters

<i>callbacks</i>	Map of callbacks per iterator
<i>block</i>	Shall this method block until all callbacks are finished

7.26.3.30 setConditionalFilter()

```
virtual void TimeTagger::setConditionalFilter (
    std::vector< channel_t > trigger,
    std::vector< channel_t > filtered ) [pure virtual]
```

configures the conditional filter

After each event on the trigger channels, one event per filtered channel will pass afterwards. This is handled in a very early stage in the pipeline, so all event limitations but the deadline are suppressed. But the accuracy of the order of those events is low.

Refer the Manual for a description of this function.

Parameters

<i>trigger</i>	the channels that sets the condition
<i>filtered</i>	the channels that are filtered by the condition

7.26.3.31 setDeadtime()

```
virtual timestamp_t TimeTagger::setDeadtime (
    channel_t channel,
    timestamp_t deadtime ) [pure virtual]
```

set the deadtime between two edges on the same channel.

This function sets the user configurable deadtime. The requested time will be rounded to the nearest multiple of the clock time. The deadtime will also be clamped to device specific limitations.

As the actual deadtime will be altered, the real value will be returned.

Parameters

<i>channel</i>	channel to be configured
<i>deadtime</i>	new deadtime

Returns

the real configured deadtime

7.26.3.32 setEventDivider()

```
virtual void TimeTagger::setEventDivider (
    channel_t channel,
    unsigned int divider ) [pure virtual]
```

Divides the amount of transmitted edge per channel.

This filter decimates the events on a given channel by a specified. factor. So for a divider n, every nth event is transmitted through the filter and n-1 events are skipped between consecutive transmitted events. If a conditional filter is also active, the event divider is applied after the conditional filter, so the conditional is applied to the complete event stream and only events which pass the conditional filter are forwarded to the divider.

As it is a hardware filter, it reduces the required USB bandwidth and CPU processing power, but it cannot be configured for virtual channels.

Parameters

<i>channel</i>	channel to be configured
<i>divider</i>	new divider, must be smaller than 65536

7.26.3.33 setFilter()

```
virtual void TimeTagger::setFilter (
    bool state ) [pure virtual]
```

enables or disables the filter on the FPGA board.

Deprecated use setConditionalFilter

The filter disables transmission of nontriggered tags on channel 7. This is a deprecated specialization of setConditionalFilter.

7.26.3.34 setHardwareBufferSize()

```
virtual void TimeTagger::setHardwareBufferSize (
    int size ) [pure virtual]
```

sets the maximum USB buffer size

This option controls the maximum buffer size of the USB connection. This can be used to balance low input latency vs high (peak) throughput.

Parameters

<i>size</i>	the maximum buffer size in events
-------------	-----------------------------------

7.26.3.35 setInputDelay()

```
virtual void TimeTagger::setInputDelay (
    channel_t channel,
    timestamp_t delay ) [pure virtual]
```

set time delay on a channel

When set, every event on the channel is delayed by the given delay in picoseconds. Setting larger time delays consumes, dependend on input signal, a significant amount of memory. Implementers are adviced to keep the delay below 1 micro second.

Parameters

<i>channel</i>	the channel to set
<i>delay</i>	the delay in picoseconds

7.26.3.36 setLED()

```
virtual void TimeTagger::setLED (
    uint32_t bitmask ) [pure virtual]
```

Enforce a state to the LEDs 0: led_status[R] 16: led_status[R] - mux 1: led_status[B] 17: led_status[B] - mux 2: led_status[G] 18: led_status[G] - mux 3: led_power[R] 19: led_power[R] - mux 4: led_power[B] 20: led_power[B] - mux 5: led_power[G] 21: led_power[G] - mux 6: led_clock[R] 22: led_clock[R] - mux 7: led_clock[B] 23: led_clock[B] - mux 8: led_clock[G] 24: led_clock[G] - mux.

7.26.3.37 setNormalization()

```
virtual void TimeTagger::setNormalization (
    bool state ) [pure virtual]
```

enables or disables the normalization of the distribution.

Refer the Manual for a description of this function.

7.26.3.38 setTestSignal()

```
virtual void TimeTagger::setTestSignal (
    channel_t channel,
    bool enabled ) [pure virtual]
```

enable the calibration on a channel.

This will connect or disconnect the channel with the on-chip uncorrelated signal generator.

Parameters

<i>channel</i>	the channel
<i>enabled</i>	enabled / disabled flag

7.26.3.39 setTestSignalDivider()

```
virtual void TimeTagger::setTestSignalDivider (
    int divider ) [pure virtual]
```

set the divider for the frequency of the test signal

The base clock of the test signal oscillator for the Time Tagger Ultra is running at 100.8 MHz sampled down by an factor of 2 to have a similar base clock as the Time Tagger 20 (~50 MHz). The default divider is 63 -> ~800 kEvents/s

Parameters

<i>divider</i>	frequency divisor of the oscillator
----------------	-------------------------------------

7.26.3.40 setTriggerLevel()

```
virtual void TimeTagger::setTriggerLevel (
    channel_t channel,
    double voltage ) [pure virtual]
```

set the trigger voltage threshold of a channel

Parameters

<i>channel</i>	the channel to set
<i>voltage</i>	voltage level.. [0..1]

7.26.3.41 sync()

```
virtual void TimeTagger::sync ( ) [pure virtual]
```

Sync the timetagger pipeline, so that all started iterators and their enabled channels are ready.

7.26.3.42 unregisterChannel()

```
virtual void TimeTagger::unregisterChannel (
    channel_t channel ) [pure virtual]
```

release a previously registered channel.

Parameters

<i>channel</i>	the channel
----------------	-------------

The documentation for this class was generated from the following file:

- TimeTagger.h

7.27 TimetaggerFPGA Class Reference

Public Member Functions

- **TimetaggerFPGA** (std::string serial="", int blocksize=1024)
- int **configure** ()
- bool **configured** ()
- void **setTestSignalDivider** (int divider)
- void **setTriggerLevel** (int channel, double voltage)
- bool **sendDacCommand** (int prefix, int control, int address, int data, int feature)
- int **read** (void *buffer, int buffersize)
- bool **open** (std::string serial)
- std::string **getSerial** ()
- std::string **getModel** ()
- std::vector< double > **getDACRange** ()
- bool **setWireIn** (unsigned addr, int value)
- int **getWireOut** (unsigned addr)
- long **WriteToPipeIn** (int epAddr, long length, void *data)
- bool **ActivateTrigger** (unsigned addr, unsigned bit)
- bool **isTriggered** (unsigned addr, unsigned bit)
- bool **UpdateWireIns** ()
- bool **UpdateWireOuts** ()
- bool **UpdateTriggerOuts** ()
- std::string **getSensorData** ()
- std::vector< char > **SendToFlash** (const char *command, int input_length, int output_length=0, bool blocking=false)
- void **setLED** (uint32_t bitmask)

Static Public Member Functions

- static std::vector< std::string > **getDeviceList** ()
- static std::string **getTimeTaggerModel** (std::string serial)

The documentation for this class was generated from the following files:

- TimetaggerFPGA.h
- TimetaggerFPGA.cpp

7.28 TimeTaggerModel Class Reference

Static Public Member Functions

- static std::string **mapOkModelToTimeTaggerModel** (OpalKellyLegacy::okCFrontPanel::BoardModel model)

Static Public Attributes

- static const std::string **MODEL_UNKNOWN** = "unknown"
- static const std::string **MODEL_TIMETAGGER_20** = "Time Tagger 20"
- static const std::string **MODEL_TIMETAGGER_ULTRA** = "Time Tagger Ultra"

Friends

- class **TimetaggerFPGA**

The documentation for this class was generated from the following files:


- TimetaggerFPGA.h
- TimetaggerFPGA.cpp

7.29 TimeTagStream Class Reference

access the time tag stream

```
#include <Iterators.h>
```

Inheritance diagram for TimeTagStream:



Public Member Functions

- [TimeTagStream](#) ([TimeTagger](#) *tagger, int n_max_events, const std::vector< channel_t > &channels=std::vector< channel_t >())
constructor of a [TimeTagStream](#) thread
- [~TimeTagStream](#) ()
tbd
- size_t [getCounts](#) ()
get incoming time tags
- void [getData](#) ([TimeTagStreamBuffer](#) &buffer)

Protected Member Functions

- bool [next_impl](#) (std::vector< [Tag](#) > &incoming_tags, timestamp_t begin_time, timestamp_t end_time) override
update iterator state
- void [clear_impl](#) () override
clear [Iterator](#) state.

Additional Inherited Members

7.29.1 Detailed Description

access the time tag stream

7.29.2 Constructor & Destructor Documentation

7.29.2.1 ~TimeTagStream()

```
TimeTagStream::~TimeTagStream ( )
```

tbd

7.29.3 Member Function Documentation

7.29.3.1 clear_impl()

```
void TimeTagStream::clear_impl ( ) [override], [protected], [virtual]
```

clear [Iterator](#) state.

Each [Iterator](#) should implement the [clear_impl\(\)](#) method to reset its internal state. The [clear_impl\(\)](#) function is guarded by the update lock.

Reimplemented from [_Iterator](#).

7.29.3.2 getCounts()

```
size_t TimeTagStream::getCounts ( )
```

get incoming time tags

All incoming time tags are stored in a buffer (max size: max_tags). The buffer is cleared after retrieving the data with getData() return the number of stored tags

7.29.3.3 next_impl()

```
bool TimeTagStream::next_impl (
    std::vector< Tag > & incoming_tags,
    timestamp_t begin_time,
    timestamp_t end_time ) [override], [protected], [virtual]
```

update iterator state

Each [Iterator](#) must implement the [next_impl\(\)](#) method. The [next_impl\(\)](#) function is guarded by the update lock.

The backend delivers each [Tag](#) on each registered channel to this callback function.

Parameters

<i>list</i>	block of events
<i>start_time</i>	earliest event in the block
<i>end_time</i>	start_time of the next block, not including in this block

Returns

if the content of this block was modified

Implements [_Iterator](#).

The documentation for this class was generated from the following files:

- Iterators.h
- Iterator.cpp

7.30 TimeTagStreamBuffer Class Reference

Public Member Functions

- **GET_DATA_1D** (getOverflows, unsigned char, array_out,)
- **GET_DATA_1D** (getChannels, channel_t, array_out,)
- **GET_DATA_1D** (getTimestamps, timestamp_t, array_out,)

Public Attributes

- `std::vector< unsigned char >` **tagOverflows**
- `std::vector< channel_t >` **tagChannels**
- `std::vector< timestamp_t >` **tagTimestamps**
- `int` **size**
- `bool` **hasOverflows**
- `timestamp_t` **tStart**
- `timestamp_t` **tGetData**

The documentation for this class was generated from the following file:

- `Iterators.h`

Index

- [_Iterator](#), [15](#)
 - [_Iterator](#), [17](#)
 - [~_Iterator](#), [17](#)
 - [channels_registered](#), [20](#)
 - [clear](#), [17](#)
 - [clear_impl](#), [17](#)
 - [getCaptureDuration](#), [17](#)
 - [getNewVirtualChannel](#), [18](#)
 - [isRunning](#), [18](#)
 - [lock](#), [18](#)
 - [next_impl](#), [18](#)
 - [registerChannel](#), [19](#)
 - [running](#), [20](#)
 - [start](#), [19](#)
 - [startFor](#), [19](#)
 - [stop](#), [20](#)
 - [unlock](#), [20](#)
 - [unregisterChannel](#), [20](#)
- [~Dump](#)
 - [Dump](#), [43](#)
- [~TimeTagStream](#)
 - [TimeTagStream](#), [83](#)
- [~TimeTagger](#)
 - [TimeTagger](#), [70](#)
- [~_Iterator](#)
 - [_Iterator](#), [17](#)
- [Array< T >](#), [21](#)
- [autoCalibration](#)
 - [TimeTagger](#), [70](#)
- [base iterators](#), [12](#)
 - [Dump](#), [13](#)
 - [TimeTagStream](#), [13](#)
- [channel](#)
 - [Tag](#), [63](#)
- [channels_registered](#)
 - [_Iterator](#), [20](#)
- [clear](#)
 - [_Iterator](#), [17](#)
 - [SynchronizedMeasurements](#), [62](#)
- [clear_impl](#)
 - [_Iterator](#), [17](#)
 - [Combiner](#), [26](#)
 - [Correlation](#), [29](#)
 - [CountBetweenMarkers](#), [32](#)
 - [Counter](#), [36](#)
 - [Countrate](#), [38](#)
 - [Dump](#), [43](#)
 - [HistogramLogBins](#), [52](#)
 - [Iterator](#), [54](#)
 - [Scope](#), [57](#)
 - [StartStop](#), [60](#)
 - [TimeDifferences](#), [65](#)
 - [TimeTagStream](#), [83](#)
- [clearOverflows](#)
 - [TimeTagger](#), [71](#)
- [Coincidence](#), [21](#)
 - [Coincidence](#), [22](#)
 - [getChannel](#), [22](#)
- [Coincidences](#), [23](#)
 - [Coincidences](#), [24](#)
 - [getChannels](#), [24](#)
 - [next_impl](#), [24](#)
- [Combiner](#), [25](#)
 - [clear_impl](#), [26](#)
 - [Combiner](#), [26](#)
 - [GET_DATA_1D](#), [26](#)
 - [getChannel](#), [26](#)
 - [next_impl](#), [27](#)
- [Correlation](#), [27](#)
 - [clear_impl](#), [29](#)
 - [Correlation](#), [28](#)
 - [GET_DATA_1D](#), [29](#)
 - [next_impl](#), [30](#)
- [CountBetweenMarkers](#), [30](#)
 - [clear_impl](#), [32](#)
 - [CountBetweenMarkers](#), [31](#)
 - [GET_DATA_1D](#), [32, 33](#)
 - [next_impl](#), [33](#)
 - [ready](#), [33](#)
- [Counter](#), [34](#)
 - [clear_impl](#), [36](#)
 - [Counter](#), [34](#)
 - [GET_DATA_2D](#), [36](#)
 - [next_impl](#), [36](#)
 - [start](#), [37](#)
- [Countrate](#), [37](#)
 - [clear_impl](#), [38](#)
 - [Countrate](#), [38](#)
 - [GET_DATA_1D](#), [38](#)
 - [next_impl](#), [39](#)
 - [start](#), [39](#)
- [DelayedChannel](#), [40](#)
 - [DelayedChannel](#), [40](#)
 - [getChannel](#), [41](#)
 - [next_impl](#), [41](#)
 - [setDelay](#), [41](#)

- start, 42
- Dump, 42
 - ~Dump, 43
 - base iterators, 13
 - clear_impl, 43
 - next_impl, 43
- Event, 44
- Flim, 44
 - Flim, 45
- FrequencyMultiplier, 45
 - FrequencyMultiplier, 46
 - next_impl, 47
- GET_DATA_1D
 - Combiner, 26
 - Correlation, 29
 - CountBetweenMarkers, 32, 33
 - Countrate, 38
- GET_DATA_2D
 - Counter, 36
 - TimeDifferences, 66
 - TimeTagger, 71
- GatedChannel, 47
 - GatedChannel, 48
 - getChannel, 49
 - next_impl, 49
- getCaptureDuration
 - _iterator, 17
- getChannel
 - Coincidence, 22
 - Combiner, 26
 - DelayedChannel, 41
 - GatedChannel, 49
- getChannelNumberScheme
 - TimeTagger, 71
- getChannels
 - Coincidences, 24
 - TimeTagger, 71
- getConditionalFilterFiltered
 - TimeTagger, 72
- getConditionalFilterTrigger
 - TimeTagger, 72
- getCounts
 - TimeDifferences, 66
 - TimeTagStream, 83
- getDACRange
 - TimeTagger, 72
- getDeadtime
 - TimeTagger, 72
- getEventDivider
 - TimeTagger, 73
- getFilter
 - TimeTagger, 73
- getHardwareBufferSize
 - TimeTagger, 73
- getInputDelay
 - TimeTagger, 73
- getInvertedChannel
 - TimeTagger, 74
- getModel
 - TimeTagger, 74
- getNewVirtualChannel
 - _iterator, 18
- getNormalization
 - TimeTagger, 74
- getOverflows
 - TimeTagger, 74
- getOverflowsAndClear
 - TimeTagger, 74
- getPcbVersion
 - TimeTagger, 74
- getPsPerClock
 - TimeTagger, 75
- getSensorData
 - TimeTagger, 75
- getSerial
 - TimeTagger, 75
- getTestSignal
 - TimeTagger, 75
- getTriggerLevel
 - TimeTagger, 75
- Histogram, 49
 - Histogram, 50
- HistogramLogBins, 51
 - clear_impl, 52
 - HistogramLogBins, 51
 - next_impl, 52
- isRunning
 - _iterator, 18
- isUnusedChannel
 - TimeTagger, 76
- Iterator, 53
 - clear_impl, 54
 - Iterator, 54
 - next, 54
 - next_impl, 54
 - size, 55
- LED, 55
- lock
 - _iterator, 18
- next
 - Iterator, 54
- next_impl
 - _iterator, 18
 - Coincidences, 24
 - Combiner, 27
 - Correlation, 30
 - CountBetweenMarkers, 33
 - Counter, 36
 - Countrate, 39
 - DelayedChannel, 41
 - Dump, 43

- FrequencyMultiplier, 47
 - GatedChannel, 49
 - HistogramLogBins, 52
 - Iterator, 54
 - Scope, 57
 - StartStop, 60
 - TimeDifferences, 66
 - TimeTagStream, 84
- overflow
- Tag, 63
- PRBS, 56
- ready
- CountBetweenMarkers, 33
 - TimeDifferences, 67
- registerChannel
- _Iterator, 19
 - TimeTagger, 76
- registerMeasurement
- SynchronizedMeasurements, 62
- reset
- TimeTagger, 76
- runSynchronized
- TimeTagger, 76
- running
- _Iterator, 20
- Scope, 56
- clear_impl, 57
 - next_impl, 57
 - Scope, 57
 - start, 58
- setConditionalFilter
- TimeTagger, 77
- setDeadtime
- TimeTagger, 77
- setDelay
- DelayedChannel, 41
- setEventDivider
- TimeTagger, 78
- setFilter
- TimeTagger, 78
- setHardwareBufferSize
- TimeTagger, 78
- setInputDelay
- TimeTagger, 79
- setLED
- TimeTagger, 79
- setMaxCounts
- TimeDifferences, 67
- setNormalization
- TimeTagger, 79
- setTestSignal
- TimeTagger, 79
- setTestSignalDivider
- TimeTagger, 80
- setTriggerLevel
- TimeTagger, 80
- size
- Iterator, 55
- start
- _Iterator, 19
 - Counter, 37
 - Countrate, 39
 - DelayedChannel, 42
 - Scope, 58
 - StartStop, 60
 - SynchronizedMeasurements, 62
 - TimeDifferences, 67
- startFor
- _Iterator, 19
- StartStop, 58
- clear_impl, 60
 - next_impl, 60
 - start, 60
 - StartStop, 59
- stop
- _Iterator, 20
 - SynchronizedMeasurements, 62
- sync
- TimeTagger, 80
- SynchronizedMeasurements, 61
- clear, 62
 - registerMeasurement, 62
 - start, 62
 - stop, 62
 - SynchronizedMeasurements, 61
- Tag, 63
- channel, 63
 - overflow, 63
 - time, 63
- time
- Tag, 63
- TimeDifferences, 64
- clear_impl, 65
 - GET_DATA_2D, 66
 - getCounts, 66
 - next_impl, 66
 - ready, 67
 - setMaxCounts, 67
 - start, 67
 - TimeDifferences, 65
- TimeTagStream, 82
- ~TimeTagStream, 83
 - base iterators, 13
 - clear_impl, 83
 - getCounts, 83
 - next_impl, 84
- TimeTagStreamBuffer, 84
- TimeTagger, 68
- ~TimeTagger, 70
 - autoCalibration, 70
 - clearOverflows, 71
 - GET_DATA_2D, 71
 - getChannelNumberScheme, 71

- getChannels, [71](#)
- getConditionalFilterFiltered, [72](#)
- getConditionalFilterTrigger, [72](#)
- getDACRange, [72](#)
- getDedtime, [72](#)
- getEventDivider, [73](#)
- getFilter, [73](#)
- getHardwareBufferSize, [73](#)
- getInputDelay, [73](#)
- getInvertedChannel, [74](#)
- getModel, [74](#)
- getNormalization, [74](#)
- getOverflows, [74](#)
- getOverflowsAndClear, [74](#)
- getPcbVersion, [74](#)
- getPsPerClock, [75](#)
- getSensorData, [75](#)
- getSerial, [75](#)
- getTestSignal, [75](#)
- getTriggerLevel, [75](#)
- isUnusedChannel, [76](#)
- registerChannel, [76](#)
- reset, [76](#)
- runSynchronized, [76](#)
- setConditionalFilter, [77](#)
- setDedtime, [77](#)
- setEventDivider, [78](#)
- setFilter, [78](#)
- setHardwareBufferSize, [78](#)
- setInputDelay, [79](#)
- setLED, [79](#)
- setNormalization, [79](#)
- setTestSignal, [79](#)
- setTestSignalDivider, [80](#)
- setTriggerLevel, [80](#)
- sync, [80](#)
- TimeTagger, [70](#)
- unregisterChannel, [81](#)
- TimeTagger backend, [11](#)
- TimeTaggerModel, [82](#)
- TimetaggerFPGA, [81](#)
- unlock
 - _Iterator, [20](#)
- unregisterChannel
 - _Iterator, [20](#)
 - TimeTagger, [81](#)